



# Daniel Mini

Paediatric bed  
with integrated weighing system



BED WITH  
INTEGRATED SCALE



**9 NOVELTY**



total  
**VISION**  
Safety and well-being



[www.gardhenbilance.it](http://www.gardhenbilance.it)



standard equipment


# 9 NOVELTY

## total VISION

The transparent plexiglass bed shoulders of the "Total Vision" system ensure continuous monitoring and a facilitated interaction between the child and health professionals, thus improving the effectiveness and safety of therapy.

Side rails' cushioned descent and assisted raise

Display - GB Touch screen movements and weighing system handling through a single panel




Children are safe and protected and not in a cage

See-through head and foot panels to keep the little patients always under control

Bed with integrated scale







# 9 NOVELTY

accessories

## LIGHT EFFECT

Head and foot panels equipped with colorful lights with different intensity. Lights brightness can be adjusted through the Touch screen



## Shoulder customization

Accessory available with the following subjects.




Customizations can also be made on request - Equipment code AC185/4.





standard equipment

# 9 NOVELTY




7

**Motorized bed extension.**  
It increases bed length up to 400mm



8

**Seat/Chair Cardiac Chair**



9

**Roto-translation of the backrest**  
to avoid chest compression

Bed with integrated scale







## Centralized brake and directional system


The central brake is activated by a pedal which is easy accessible, allowing the wheels to brake simultaneously.

It is also possible to activate the directional system to facilitate the transport of the device during movements from a place to another (Equipment code AC7L2G).



### Wheels

No.4 stainless, swivable and non-marking wheels. Diam 150mm.



### Automatic cable reel

Equipment code AC50L/1.



### Hooks for bags

Located on the frame.






# Movements

## Backrest quick release

The bed is equipped with a backrest's manual quick release function (AC8L). To use this function the user will pull one of the 2 levers located at the head end. The levers are two in order to permit the user to perform this important operation from both sides of the bed. In addition the CPR position can be also achieved automatically with a single touch through the GB Touch digital display.



### CPR Control

Manual activation through lever

### Automatic CPR control

The command can be activated from the display - GBTouch Console.






# Movements

## Motorized bed extension

The bed is provided with a motorized extension that increases the overall length of the device up to 400mm more. This feature is intended to allow the bed to be suitable also for patients of bigger size. This movement is controlled through the GB Touch digital display (Equipment Code AC61/4).


### Motorized bed extension

Electrically activated through Touch screen display - GB Touch



### Additional mattresses

The equipment includes n. 2 additional 200 mm mattresses.





# Movements



## Trendelenburg and Reverse Trendelenburg

Trendelenburg and Reverse Trendelenburg are achievable by pushing on the dedicated push-buttons located on the display - GB Touch screen.

Trendelenburg



Reverse Trendelenburg








# Movements

## Adjustable height

The paediatric bed is height adjustable. The height adjustment is reached through no.2 rectangular telescopic motorized columns that are activated through the GB Touch display. The height can be adjusted up to +400mm from lowest position.






# Movements

## Cardiac chair

The bed can achieve the "cardiac chair" position. Also in this case the backrest section performs automatic roto-translation. The function is memorized and controlled via the GB Touch display console. Equipment code AC62.



**Cardiac chair position**





## GB Touch


### Scale display and Control box

The pediatric bed is equipped with a 10.1" HD Touch Screen Display with a 160° visual angle to manage the **integrated weighing system** and all the **movements** of the device.



#### Scale display and Control box

equipped with a graphical interface to manage the weighing system.



In addition, the display has been designed to be embedded into the bed's frame in order to be always visible.





## Features

- Switch-on via Soft Touch Button
- Movements management
- Low voltage power supply
- Indicator box in plastic material with IP65 insulation
- Battery status indicator
- Display lighting control
- Screen saver
- Night mode
- Smart functions
- Complete control via GBSmart android/ios/windows application
- Scheduled maintenance setting
- Language setting: Italian - English - French - Spanish - German
- Digital protection against accidental or unwanted access
- Shoulder lighting management (only when accessory AC185 has been equipped on the bed)



## Features and functionality of the weighing system

- Weight reading via 10.1" HD touch screen display
- Patient's biographical data
- Patient's location indication
- Tare: automatic and manual. This value is always displayed on the screen
- Unit of measurement indicator: unit of measurement can be selected between kg and lb
- Capacity 200kg
- Minimum division double scale automatic display:
  - 50g up to 150 kg and 100g from 150 kg up to 200 kg.
- Converter: 24bit, 10,000 divisions.
- Weight loss alarms.
- Patient descent alarms.
- Access register.
- Remote parameter control
- Automatic data transmission to enabled ports
- Ethernet port 10/100 - Equipment code AC167/3.
- Upon request: Combo WI-Fi/Bluetooth module, Accessory code AC167/1 - USB A port for weight data transmission, Accessory code AC167/4 - USB B port for external printer connection, Accessory code AC167/5 - RS-232 port - RS-485 port, Accessory code AC167/6. Battery buffer to ensure autonomy of the weighing system in the absence of power - equipment code AC17/L.



### For more information about the GB Touch

Frame or click the QR code button and discover all the features and capabilities of the display.



## Other Features

- Frame painted with thermosetting epoxy powders
- Bed base made of 3-section in HPL with double articulation
- Thermoformed ABS cover to protect the base
- Easily removable transparent shoulders with locking/unlocking mechanism
- Transparent side rails with cushioned descent and assisted raise, with release lever and integrated handle for movement
- Motorized backrest section
- Motorized femoral section
- Servo-assisted leg section by gas springs, synchronized with the movement of the femoral section and manually adjustable
- Adjustable height, motorized Trendelenburg and reverse Trendelenburg through electric columns
- Cardiological position with roto-translation of the backrest
- Motorized bed extender with up to +400 mm
- Anti-crushing and anti-shear low voltage electric motors
- Swivel, stainless, non-marking wheels, diameter 150 mm with centralized braking system and directional system. Activated by foot pedal
- Quick release of the backrest for easy CPR
- Battery box - to power the device in case of power failure
- No.4 bumpers.
- Automatic cable reel
- Equipotential node
- Predisposition for IV pole
- Bilateral bag holder hooks
- Clearance allows lifts or hoists fitting
- Touch screen display for managing movements and the weighing system - Equipment code AC167.
- Additional 10-key push-button panel. The supplied hand control is only intended for use in the event of a failure of the GB Touch handling system.
- The control panel is equipped with a key for inhibition of controls (equipment code AC34L).

### Do you want to choose the colour of the frame?


Frame or click the QR code and find out how to customize the device by consulting our color cards\*.

\*If no preference is expressed, the device will be supplied with the available color.

Frame



## Dimensions



## Main cord



**Type F**  
Schuko



**Type L**  
16A



**Type I**  
Australia



**Type B**  
Canada and USA



**Type G**  
UK plug





## Technical sheet


|                             |  |
|-----------------------------|--|
| CND - EMDN code             | <b>V0806010702</b>   |
| GMDN code                   | <b>35321</b>   |
| N. progressivo/R            | <b>2525609/R</b>   |
| UDI-DI di base              | <b>805771740LETBGV</b>   |
| ID Product                  | <b>GB0021.MP</b> Painted steel   |
| Inteded use                 | Bed with integrated scale for weight determination of the patient for assessment, diagnosis and care purposes.   |
| Manufacturer                | GARDHEN BILANCE SRL  |
| Spare parts availability    | 10 years   |
| Overall weight              | 168 kg   |
| Controls                    | Display touch screen   |
| Battery                     | 10 cycles of autonomy  |
| Battery weighing system     | 24 V self rechargeable (autonomy approx. 24 hours)   |
| Power                       | Multivoltage 100-240 Vac 50/60 Hz  |
| Main cable                  | Schuko or other upon request   |
| Wheels                      | No. 4 swiveling and stainless 150 mm with braking system   |
| Electric motors             | Low voltage (24 V)   |
| Max electrical input Max    | 460 VA   |
| Insulation                  | IP44   |
| Max load                    | 200 kg   |
| Safety working load         | 235 kg   |
| Divisions                   | With a max load capacity of 200 kg automatic displaying of the minimal division of 50g up to 150 kg and 100g from 150 kg up to 200 kg..  |
| Weighing cell power supply: | 5Vcc (from 1 to 4)   |
| Display box                 | Display box in plastic material with IP65 insulation   |
| Converter                   | 24bit  |
| Display                     | Touch Screen   |
| Power display               | 24 Vdc   |
| Alarms                      | Status indicator Luminous on display and on control systems<br>Alarms Excessive or insufficient weight loss / Patient exit alarm /<br>Acoustic brake disconnection alarm / Battery status indicator /<br>Faults indicator / Maintenance deadline alarm |
| Data transmission           | Ethernet port 10/100<br>On request : WI-Fi/Bluetooth combo module - USB A port for weight-given transmission - USB B port for external printer connection - RS-232 port - RS-485 port.   |
| I/O (opz.)                  | 4 output   |
| Antiblack-out function      | Extra battery + software   |
| Tare                        | Tare subtraction before and during the treatment   |
| Tare modalities             | Automatic and manual   |
| Conformity                  | Regulation UE 2017/745<br>Certificate No. 126/MDR issued by O.B. 0051<br>It meets all applicable standards and directives  |



# Daniel Mini

Paediatric bed with integrated weighing system

[www.gardhenbalance.it](http://www.gardhenbalance.it)



**LIGHT  
EFFECT**



## Gardhen Balance S.r.l.

Via G. Luraghi c/o Consorzio Il Sole - lotto S

80038 Pomigliano d'Arco (NA) - Italy

Tel. +39.081.8692160 pbx - Fax +39.081.8692460

E-mail: [info@gardhenbalance.it](mailto:info@gardhenbalance.it)

PEC: [gardhenbalance@pec.it](mailto:gardhenbalance@pec.it)



Pictures contained in this brochure is for illustrative purpose only. The manufacturer reserves the right to modify the projects without notice. For accessories identification please refer to the accessory catalogue or commercial offer. When data are different please refer to commercial offer.

Dep. EN 110 rev. 01/10/2024