

Straight Up Esense

PASSIVE Lift



ASSISTED MOBILITY

Frequent use

Home care
Hospital care
Retirement house
Nursing home
Rehabilitation
Long stay



LIFT UP



**SOLID &
EASY TO HANDLE**

for lifting people who are not able
to support their own weight



gardhen
balance

www.gardhenbalance.it



Standard equipment


The electric lifter **Straight Up - PASSIVE** it is suitable for lifting people who **are unable to support their own weight** and to ensure safety in **movement from a device to another** (bed, wheelchair, stretcher...) avoiding any effort.

Ergo-feedback
to prevent
operator's
physical
overloads

Handlebar
equipped with
pushing sensor

Twin wheels
of which two with brakes to
facilitate movements both
in hospitals
and at home

E-drive system
allows you to move the lift
effortlessly in any
direction even with
a patient on board





Standard equipment



125 mm front wheels with brake



Handset
7-key handset for operating the motors



Charger



Rechargeable battery



Emergency button
Emergency button to lock movements



Manual adjustment
Manual emergency lever for lifting adjustment





The electric lifter **Straight Up - PASSIVE** is designed to allow the operator to independently manage the daily activities of **transfer and repositioning the patients** from a rooms to another.

Thanks to the principle of vertical lifting and the **width adjustment of the legs of the base**, it is ideal for use in domestic and residential environments and in the rooms where the **space under the bed is limited**.

The set of twin wheels with brake that can be activated by foot pedal, the practical activation control with up/down button allow you to **move patients who would be difficult to move manually**.





Ergo-feedback

The electric lifter **Straight Up - PASSIVE** is equipped with an Ergo-feedback system. This safety system helps operators to prevent physical overloads as it advises in real time if the lift is moving with too much force (**Orange LEDs that flash**)

Flashing warning thrust force control



Activation and Orientation Automatic

E-drive

The electric lifter **Straight Up - PASSIVE PASSIVO** is equipped with an E-drive system which allows you to move the lift effortlessly in any direction even with the patient on board.

The handlebar is equipped with a sensor which detects the thrust and communicates the direction to the 5th wheel.

The presence of the E-drive system in the lift facilitates the work of the operators, avoiding efforts that in the long run could cause body pain.

The activation of the E-drive system occurs by automatically orienting the wheel.

Easily maneuverable by applying a light push





Movements - Lifting

The **lifting and lowering** of the patient is controlled through to dedicated keys on the handset.



Lift Up



Low Down





Movements - Inclination of the slingbar for the harness

The lifter is equipped with sling bar and a electric 4 point spreaderbar. The **inclination of the sling bar** is controlled trough the dedicated buttons on the **handset**, allowing to adjust the patient from lying to sitting.



Lying down




Sitting





Movements - Legs of the device

The lifter has a mechanism that allows you to adjust the **width of the legs of the base of the device** to bring it closer to the patient if sitting on a chair, armchair or lying on a bed. The width of the legs can be full **closed** for moving or positioning under the bed. Movement is controlled trough the dedicated keys on the **handset**.



closed legs for movement or for positioning under the bed

opening of the legs
to get closer
to the patient
when seated



Opening



Closure





Standard equipment

- Steel frame painted with thermosetting powders.
- Operator handlebar with pushing sensor
- Nr. 2 front twin wheels with brake activated by pedal - 125 mm
- Nr. 2 twin wheels which allow easy accessibility under the bed or other devices - 100 mm
- Handset equipped with coiled cable. The handset is easy to use as it is made with large relief pictograms and is also suitable for visually impaired patients.
- E-drive system
- Ergo-feedback
- Manual emergency lever for lifting adjustment
- Rechargeable battery
- Software
- Output connector for data analysis
- Capacity 175 kg.

Data sheet

Product code	GBSL0001 - Straight Up Esense GBSL0002 - Straight Up
Intended use	Used to lift patients who are unable to support their own weight

Power cable type



Type L - 16A



Accessories

Slings

Lifting slings are used in numerous care activities including lifting, positioning and transferring patients from beds, chairs, etc.

Available in different sizes.

SEAMLESS Sling

Model that **guarantees support for the whole body**. It features a comfortable headsupport to prevent inappropriate cervical movements. The sling has no seams and prevents pressure ulcers. In order to position this sling, the patient must either be lying in bed



STANDARD Sling

Model available in 2 fabrics, **black polyester** and **black net polyester** for use in swimming pools and bathrooms.

The sling can be used when the patient is sitting in the chair but not in the sling.

The sling is positioned up to the patient's sacrum and not under the buttocks. When the patient is in the chair, the sling can be removed



	SUPPORT WHOLE BODY	WEARABLE DA WHEN SITTING	WITHOUT SEAMS	SUPPORT FOR HEADSUPPORT	TRANSFER IN THE POOL AND BATH
 SEAMLESS Sling	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
 STANDARD Sling	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>



STANDARD Sling - Net





The sling attachment system consists of **fixing clips** which signal the locking or unlocking of the sling to the lifter.



Harness attachment clip

If the red color appears on the clip the locking has not been carried out, if instead it is of green color the clip is perfectly fixed to the lifter. The universal clip is also suitable for other lifts on the market according to the ISO-10535 standard.

With the safety lock mechanism



RED - UNLOCKED



GREEN - BLOCKED

Sling - Sizes available

Size	Seamless Sling	Standard Sling	Standard Mesh Sling
XXS	GBSL0027	GBSL0009	GBSL0018
XS	GBSL0028	GBSL0010	GBSL0019
S	GBSL0029	GBSL0011	GBSL0020
M	GBSL0030	GBSL0012	GBSL0021
L	GBSL0031	GBSL0013	GBSL0022
LL	GBSL0032	GBSL0014	GBSL0023
XL	GBSL0033	GBSL0015	GBSL0024
XXL	GBSL0034	GBSL0016	GBSL0025
XXXL	GBSL0035	GBSL0017	GBSL0026



Straight Up

PASSIVE Lift

www.gardhenbilance.it



version without
E-drive system



Gardhen Bilance S.r.l.

Via G. Luraghi c/o Consorzio Il Sole - lotto S - 80038 Pomigliano d'Arco (NA) - Italia

Tel. +39.081.8692160 pbx - Fax +39.081.8692460

E-mail: info@gardhenbilance.it - PEC: gardhenbilance@pec.it

Visit our website www.gardhenbilance.it / Follow us on

The images reproduced in this sheet are purely indicative and the dimensional values are nominal. The manufacturer reserves the right to modify the projects without notice. For accessories identification please refer to the accessory catalogue or commercial offer. When data are different please refer to commercial offer.

Dep. EN 365 rev 15/04/2024