#### Theia 2.0

#### Italian design and technology in the world







ENTRY POSITION ACHIEVABLE
THROUGH MEMORY BUTTONS
LOCATED ON THE CONTROL PANEL.
REDUCING THE HEIGHT TO ITS
MINIMUM ALLOWS AN EASY ENTRANCE
ALSO TO PATIENT WITH REDUCED
MOBILITY.



EASY AND COMFORTABLE EGRESS
FOR PATIENT WITH REDUCED
MOBILITY. IT CAN BE ACHIEVED
THROUGH THE MEMORY BUTTONS ON
THE CONTROL PANEL.

#### **Benefits**

- Available in 2 versions, with high and low backrest.
- Multifunctional surgery chair with high stability. It can be used for a wide range of surgical interventions.
- Frame with a reduced footprint. It allows an easy approach to the device.
- Extremely resistant base in stainless steel. Easy to sanitize with most commonly used detergents.
- Lying surface in expanded foam with high density.
- Seamless upholstery made out of carefully selected fabric with antibacterial, water resistant and non-toxic. Upholstery's features as follows: antibacterial, flame retardant class 1MM, antimycotic, ecologic, latex free, phtalate free and resistant.
- Minimum height (530 mm, minimum value compared to market standards) makes it easy for any patient to access even to people with reduced mobility.
- Assisted egress position helping people with reduced mobility.
- Independent movement for the best postural adjustment.
- The armrests are synchornized with backrest's movement. They can be rotated outward, overturned and removed.
- No wear parts.
- The device can be equipped with a wide range of accessories according to surgeon's needs.
- Anti-schock position.
- H24 operativity without connection to the electrical grid.
- Max load capacity 300kg





#### Standard features

- Lying surface divided in 4 sections with double articulation divided in: headrest, backrest, seat and leg.
- Motorized backrest.
- Motorized legrest.
- Height adjustble through motorized columns.
- Motorized Trendelenburg.
- Headrest available from the following: ACS1, ACS2, ACS3, ACS4, ACS4/1; pediatric models ACS1P, ACS2P, ACS3P are also available. In the absence of preferences the code ACS1 will be provided. All headrests are described and illustrated on p. 13.
- Special coating without seams or interstices, made with carefully selected technological material with the following properties: water repellent, non-toxic, antibacterial, antifungal, antistatic, ecological, latex free, without phthalates. Ultra-resistant coating to: alcohol, hydrogen peroxide, sodium hypochlorite (5%), commonly used disinfectants, liquids and physiological substances. Fireproof eco-leather UNI 9175/87 and 9175FA-1/94 class 1M. (equipment code AC87 / 1).
- Upholstery in black fabric (other colours available download colours chart via QR code at the bottom of this page). Upholstery divided into 2 separate sections, one for the back section and the other for the seat / leg sections. The operating table is also available with a single seamless cover (accessory code ACS27).
- Movements are entirely controlled by the control panel (standard equipment) and/or foot controls (accessory see page 13).
- Low voltage 4000N motors.
- Multi-voltage power supply.
- Bilateral multifunctional armrest.
- Twin wheels diameter 100 mm with independent brake.
- No.8 stainless steel DIN bars that can be equipped with a wide range of accessories.
- Non-reflecting coating with thermosetting powders.
- A dedicated software allows an easy handling of the movements.
- Equipotential node.
- No.2 rechargeable and removable batteries. They grant a h24 operativity to the device.
- Charger included.
- Knobs for transport located on the leg section.



HEADRES
ACCESSORY CODE ACS1



CONTROL PANEL



INDEPENDENT BRAKING SYSTEM



KNOBS FOR TRANSPORT LOCATED AT THE LEG END



**Would you like to customize upholstery or frame?**Scan or click the QR code and discover how to customize the device consulting our colour chart.



#### Control panel

The control panel installed on the device is intuitive and easy to use.

Thanks to the pictograms depicted on the individual buttons, the surgeon can easily identify and operate the chair functions during his / her activities.

The panel is divided into sections, each of which allows the activation of specific functions, such as:

- Section dedicated to user programmable memories.
- Section dedicated to preset memories.
- Section dedicated to movements.



#### Available memory positions

The operator can easily and intuitively program, according to need, 3 different position configurations, respectively with the buttons:

- · SET1;
- · SET2:
- · SET3.

All position configurations can be recalled by pressing the dedicated keys in "push and go" or "push and stop" mode.

It is possible to interrupt the movement of the lying surface at any time. To confirm and modify the memory positions, use the **SAVE SET** button.





**BATTERY STATUS INDICATOR** 



EXAMPLE OF A MEMORY SET POSITION



#### Memory positions "PRESET"

4 frequently used memories can be activated in the ways through the control panel. "push and go" and "push and stop".

#### Memory 1 Entry and transport position



PERFECT AS ENTRYWAY, EXIT AND FOR PATIENT'S TRANSPORT.

#### Memory 2 Egress position



COMFORTABLE ASSISTED EXIT FOR PATIENTS WITH REDUCED MOBILITY IT PROGRESSIVELY VERTICALIZES THE PATIENT.



Theia 2.0

#### Memory positions "PRESET"

#### Memory 3 Trendelenburg / Anti-Schock



IT QUICKLY ACHIEVES THE PERFECT
POSITION TO FAVOUR
VENOUS RETURN

#### Memory 4 CPR



CARDIO-PULMONARY RESUSCITATION POSITION









#### **Motorized movements**

#### **Backrest adjustment**



**BUTTON FOR BACKREST ADJUSTMENT** 

#### Legrest adjustment



**BUTTON FOR LEGREST ADJUSTMENT** 





#### **Motorized movements**

#### Seat height adjustment



BUTTON FOR LYING SURFACE'S HEIGHT REGULATION

#### Seat's inclination adjustment



BUTTON FOR SEAT'S INCLINATION ADJUSTMENT





#### Headrest movement and adjustment • • •


Multi-adjustable headrest section for surgery.

The head support can be set between different angles for a correct position of the head during the operation.

It has a unique locking system, which locks both pins using a cam lever.

The different types of headrests are interchangeable thanks to a quick coupling without the use of any tools (see accessories catalog).







THANKS TO THE LEVER LOCATED IN THE MIDDLE OF THE "X" MECHANISM, IT IS POSSIBLE TO ADJUST THE INCLINATION AND THE ANGLE OF THE HEADREST IN AN EASY AND STRAIGHTFORWARD WAY







#### **Armrest**

Anatomical armrests in PU rubber with metal core.

They follow the movement of the backrest. They can be folded down to facilitate patient entry. Completely removable to facilitate patient transfer to another device. The external regulation facilitates vascular access.



TILTABLE ARMREST TO FAVOUR PATIENT'S ENTRY FROM THE SIDE



**OUTWARD REGULATION TO FAVOUR VASCULAR ACCESS** 



REMOVABLE ARMREST



#### Continuous operational autonomy

The surgical chair is equipped with two rechargeable and interchangeable batteries, with charger included. As the batteries can be changed, permanent autonomy is fully ensured. A luminous and acoustic warning will light up in case of low battery condition. The device can be powered also through the electric grid.





PLUGGED BATTERY






#### A Pedal controls


The foot controls located on the floor allow to adjust the height and the inclination of the seat section, the back section and the leg section, according to the needs of the surgeon and the type of operations to be performed, . The floor control is also available in the Wireless version (accessory code ACS5 / 3).

#### ACS5/1




Pedal on the floor for height amd inclination adjustment of the seat. Couple.

#### ACS5/2




Pedal on the floor for height adjustment of the seat and adjust the inclination of the backrest. Couple.

#### ACS6/1



Pedals on the floor, to adjust the inclination of the backrest and the inclination of the legrest. Couple.

#### ACS5/3



Pedal on the floor with wireless controls for height amd inclination of the backrest. Couple.

#### Headrest

Given that one of the following headrests is supplied with the Theia 2.0 device (if none are preferred, the ACS1 code will be provided), other models identified with the codes ACS1, ACS2, ACS3, ACS4, ACS4 / 1 are also available. pediatric models are also available ACS1P, ACS2P, ACS3P. The headrests are all interchangeable thanks to a guick coupling system (without the use of tools).

# Occinital headrest Designed to support the






Nuchal headrest. Available also in paediatric version ACS2P.



"X" shaped headrest. Available also in peadiatric version ACS3P.




Head section on DIN bar with face hole equipped with closing cap with angle adjustment of 0°-15°.



Synchronous head section with face hole, bilateral regulation from -15° to +15°. Available with head strap (ACS37) and chin strap (ACS38).





Bilateral side rails on DIN bar. Lightweight side rails, washable, adjustable and collapsible. Easy to remove.

#### ACS8/2



Support on DIN bar for humerus surgery with stand support. To be used with patient in prone position. To be fitted onto backrest of the device for humerus surgery.

#### ACS8




Hand surgery table on DIN bar. Availablewith 2 clamps to fix the support on the DIN bar. Extremely resistant and easyto use. Radiolucent. The fabric is water repellent, no toxic, antibacterial and high resistant tobiological liquid and disinfectant.


#### ACS8/1




Hand surgery table on DIN bar withstand support. Available with clampand adjustable in height in order tobe placed at any height of the tabletop. The fabric is water repellent, non-toxic, antibacterial and high resistant tobiological liquid and disinfectant.











Multifunctional armrest on DIN bar.
Adjustable in height, swivelling, reclinable, tilting, easily removable. It is intended to be an alternative to the standard armrest.



Control panel attachment on DIN bar. The pole is flexible and can be adjusted according to surgeon's needs.

#### ACS16



Flexible support for tent. Available with two clamps. To be used during surgery in the operating room in order to protect and cover the patient.

#### ACS18




Occipital pillow for headrest.

#### ACS19



Body strap on DIN bar. Fixingsystem to allow the immobilisation of the patient during surgery. Adjustable in lenght according to patient's size.



Extra padded cushion on DIN bar to enlarge the seat section. It increase the table top surface allowing the patient to comfortably lay

#### ACS24





Table on DIN bar. Available with clamp, it is designed to place the tools needed for surgical operations.

#### ACS25



Universal stainless steel DIN clamp to attach accessories with radial setting. It has a single knob that allows both fixation onto the side rail.

#### ACS26



Air flexible rod. Quick oxygen intake for ophthalmic surgery. Flexible pole that can be


One padded cushion for table top. The table top has an upholstery consisting of one piece cushion, without any interruptions.



Slim base without wheels. The device will stay on the ground.




Porte-rouleaux sur harre DIN



Oxygen tank holder on DIN bar. It can be placed on several points of the device thanks to a quick fitting on DIN bar. Up to 7 kg.

























#### **Dimensions**





#### MAIN CORD



Type F - Schuko



Type L - 16A



Type I - Australia



Type B - Canada and USA



Type G - UK plug

#### Technical sheet

GMND Code 38447

Product ID GB0210.SP-2.0 Intended use Surgery chair

Manufacturer GARDHEN BILANCE SRL

Total weight 120 kg

Type of control Control panel

Wheels Stainless swivel 4 wheels Ø 100 mm with

braking system

Power 100-240 Vac 50/60 Hz

Battery power 24v - 4.5A

Electric motors Low voltage 24 (Volt)

Max electrical input480 VAIsulationIP44Max load300 kgSafety working load335 kgBackrest angle76°Legs angle50°



Dep. EN 310 rev. 02/09/2022