

# STAZIONE DI CARICA AG7/AG300



A Phoenix Mecano Brand

Page 1/3 - stand of 02/2020

## Ricarica in tutto il mondo per le batterie

Questa stazione automatica di carica può essere utilizzata con un numero massimo di quattro batterie e rende possibile il monitoraggio di ogni singola batteria. È idonea per AG300 e AG7. L'apparecchio offre una possibilità di carica indipendentemente dalla sede del letto medico o ospedaliero. Un indicatore ottico dello stato di carica mostra l'avanzamento della ricarica della batteria.

## Soluzione di carica per il settore medico

La stazione di carica convince per il design compatto e l'alloggiamento solido e facile da pulire, caratteristiche che la rendono ideale per le applicazioni mediche.

## Vantaggi in breve

- Design compatto
- Idonea per la ricarica ottimale di batterie AG7 e AG300
- Molto solida e resistente
- Possibilità di ricarica indipendentemente dalla posizione del letto
- Facile da pulire



Note: Responsibility for DEWERT products when used for specific applications and for adherence to the appropriate guidelines, standards and laws is borne by the manufacturer of the complete system into which DEWERT products are fitted. We accept no liability for the correctness.

DewertOkin GmbH • Weststraße 1 • 32278 Kirchlengern  
Phone: +49 (0)52 23 / 9 79 -0 • Fax: +49 (0)52 23 / 7 51 82 • Email: info@dewert.de • Internet: <http://www.dewert.de>

# STAZIONE DI CARICA AG7/AG300



A Phoenix Mecano Brand

Page 2/3 - stand of 02/2020

## Standard configurations

Colore	grigio chiaro RAL 7035
Umidità relativa	30% - 75%
Temperatura ambientale	+10°C - +40°C

Note: Responsibility for DEWERT products when used for specific applications and for adherence to the appropriate guidelines, standards and laws is borne by the manufacturer of the complete system into which DEWERT products are fitted. We accept no liability for the correctness.

DewertOkin GmbH • Weststraße 1 • 32278 Kirchlengern  
Phone: +49 (0)52 23 / 9 79 -0 • Fax: +49 (0)52 23 / 7 51 82 • Email: [info@dewert.de](mailto:info@dewert.de) • Internet: <http://www.dewert.de>


# STAZIONE DI CARICA AG7/AG300



**DEWERT**

A Phoenix Mecano Brand

Page 3/3 - stand of 02/2020



Note: Responsibility for DEWERT products when used for specific applications and for adherence to the appropriate guidelines, standards and laws is borne by the manufacturer of the complete system into which DEWERT products are fitted. We accept no liability for the correctness.

DewertOkin GmbH • Weststraße 1 • 32278 Kirchlengern

Phone: +49 (0)52 23 / 9 79 -0 • Fax: +49 (0)52 23 / 7 51 82 • Email: [info@dewert.de](mailto:info@dewert.de) • Internet: <http://www.dewert.de>