

CHAIRS

ACCESSORIES

Accessory compatibility table

1	Stephen H Anatomical	27	Athena
2	Stephen H Anatomical - Mini	28	Home Chair
3	Stephen H Anatomical - Junior	29	Stephen H Ybrid
4	Stephen H Anatomical - Bariatrica	30	Stephen H Easy
6	Stephen H Flat	31	Up
7	Stephen H Flat - Junior	32	Akka
8	Stephen H Flat - Bariatrica	33	Stephen Scale H Anatomical
9	Stephen H Freedom	34	Stephen Scale H Anatomical - Mini
10	Stephen H Freedom - Junior	35	Stephen Scale H Anatomical - Junior
11	Stephen H Freedom - Bariatrica	36	Stephen Scale H Anatomical - Bariatrica
12	Stephen H Comfort 2.0 - 3.0	37	Stephen Scale H Flat
13	Stephen H Comfort 2.0 - 3.0 - Bariatrica	38	Stephen Scale H Flat - Junior
14	Stephen H Comfort 4.0	39	Stephen Scale H Flat - Bariatrica
15	Stephen H Comfort 5.0	40	Stephen Scale H Freedom
16	Stephen H Comfort 6.0	41	Stephen Scale H Freedom - Junior
17	Stephen H Comfort 7.0	42	Stephen Scale H Freedom - Bariatrica
19	Stephen H Comfort 9.0	43	Stephen Scale H Comfort 2.0 - 3.0
20	Stephen H Eco - versione 1	44	Stephen Scale H Comfort 2.0 - 3.0 - Bariatrica
21	Stephen H Eco - versione 1 - Junior	45	Athena Scale
22	Stephen H Eco - versione 2	46	CM/6R
23	Stephen H Eco - versione 2 - Junior	47	Stephen Scale
24	Stephen H Eco 2.0	48	Stephen H Infinity 1.0
25	Unica	49	Stephen
26	Athena 2.0		

EXAMPLE OF ACCESSORY CONFIGURATION

CODE
ACCESSORY

DESCRIPTION
ACCESSORY

PHOTO
PRODUCT

AC2P




Table on rod accessories for objects, equipped with height adjustment and rotation functionality.

Easily removable if necessary.

THIS ACCESSORY IS COMPATIBLE WITH ARMCHAIRS LISTED


COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	19	20	21	22	23	24	27	29
31	33	34	35	36	37	38	39	40	41	42	43	44	45	48	49								

Once the device has been configured, consult the Gardhen Balance sales offices in order to avoid any incompatibilities between the chosen accessories.

The images reproduced are purely indicative and the manufacturer has the right to make changes to the projects without notice.

AC2P



Serving table on accessory rod equipped with height and rotation adjustment. Easily removable. Dimensions 234x185 mm.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	19	20	21	22	23	24	27	29
31	33	34	35	36	37	38	39	40	41	42	43	44	45	47	48	49							

AC3P



Hook for drainage bags also available on DIN rail.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	19	20	21	22	23	24	26	27
29	30	31	33	34	35	36	37	38	39	40	41	42	43	44	45	47	48	49					

AC4P



Removable, height-adjustable **adjustable IV pole** by mechanism.
Plastic hooks - available with 2 and 4 hooks

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	19	20	21	22	23	24	27	29
30	31	33	34	35	36	37	38	39	40	41	42	43	44	45	47	48	49						

AC4P/1



Removable, height-adjustable **non adjustable IV pole** by mechanism.
Plastic hooks - available with 2 and 4 hooks.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	19	20	21	22	23	24	27	29
31	33	34	35	36	37	38	39	40	41	42	43	44	45	47	48	49							

AC4P/2



Height-adjustable **IV stand on wheels** with 4 hooks.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	19	20	21	22	23	24	25	26
27	28	29	30	31	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49		

AC4P/3



IV pole on DIN rail.

COMPATIBLE FOR THE FOLLOWING CHAIRS

6	7	8	12	13	22	26	28	37	38	39	43	44
---	---	---	----	----	----	----	----	----	----	----	----	----

AC4P/4




Shaped steel IV pole with 4 hooks..

COMPATIBLE FOR THE FOLLOWING CHAIRS

26

AC4P/5



Removable, height-adjustable **adjustable IV pole** by mechanism.
Stainless steel hooks - available with 2 and 4 hooks.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	19	20	21	22	23	24	27	29
30	31	33	34	35	36	37	38	39	40	41	42	43	44	45	47	48	49						

AC4P/6



Removable, height-adjustable **non adjustable IV pole** by mechanism.
Stainless steel hooks - available with 2 and 4 hooks.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	19	20	21	22	23	24	27	29
31	33	34	35	36	37	38	39	40	41	42	43	44	45	47	48	49							

AC4P/7




IV pole with holder, swivelling, rotatable and height-adjustable by mechanism.
Plastic hooks - available with 2 and 4 hooks.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	19	20	21	22	23	24	27	29
30	31	33	34	35	36	37	38	39	40	41	42	43	44	45	48								

AC4P/S



Bracket for IV pole.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	12	13
14	15	16	17	19	20	21	22	23	24	26	27
28	29	30	31	33	34	35	36	37	38	39	40
41	42	43	44	45	47	48	49				



Bracket for IV pole.

COMPATIBLE FOR THE FOLLOWING CHAIRS

32

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	19	20	21	22	23	24	27	29
31	33	34	35	36	37	38	39	40	41	42	43	44	45	47	48	49							

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	20	21	22	23	24	26	27	29
33	34	35	36	37	38	39	40	41	42	43	44	45	47	49									



1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	20	21	22	23	24	26	27	29
33	34	35	36	37	38	39	40	41	42	43	44	45	47	49									



1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	20	21	22	23	24	26	27	29
33	34	35	36	37	38	39	40	41	42	43	44	45	47	49									



1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	20	21	22	23	24	26	27	29
33	34	35	36	37	38	39	40	41	42	43	44	45	47	49									



1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	20	21	22	23	24	26	27	29
33	34	35	36	37	38	39	40	41	42	43	44	45	47	49									

AC6P/5



Variable height adjustment. Pedal controls located on the base. (Configurable with AC28 accessory).

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	20	21	22	23	24	26	27	29
33	34	35	36	37	38	39	40	41	42	43	44	45	47	49									


AC6P/6



Backrest adjustment. Pedal controls located on the base.

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	20	21	22	23	24	26	27	29
33	34	35	36	37	38	39	40	41	42	43	44	45	47	49									

AC6P/7




Wireless foot control.


1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	20	21	22	23	24	26	27	29
33	34	35	36	37	38	39	40	41	42	43	44	45	47	49									




1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	20	21	22	23	24	26	27	33
34	35	36	37	38	39	40	41	42	43	44	45	47	49										



1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	20	21	22	23	24	26	27	29
33	34	35	36	37	38	39	40	41	42	43	44	45	47	49									




1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	20	21	22	23	24	26	27	29
33	34	35	36	37	38	39	40	41	42	43	44	45	47	49									




Centralized brake on 4 wheels with directional system and activation via rear bilateral pedal, in galvanized steel.

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	19	26	27	29	31	33	34	35
36	37	38	39	40	41	42	43	44	45	47	48	49											



Centralized brake on 4 wheels with directional system and operation via rear bilateral pedal, in stainless steel.

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	26	27	29	31	33	34	35	36
37	38	39	40	41	42	43	44	45	47	48	49												



Bilateral brake pedal front and rear.

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	19	26	27	29	31	33	34	35
36	37	38	39	40	41	42	43	44	45	47	48	49											

AC9B



Shaped and sewn armrest cover made of 50 g/
m² TNT.

Pack of 2 pieces, minimum order quantity 100 pieces.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	19	20	21	22	23	24	25	26
27	28	29	30	31	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49		

AC9P



Shaped and sewn armchair cover made of 50 g/m2 TNT. Pack of 1 piece, minimum order quantity 60 pieces.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	19	20	21	22	23	24	25	26
27	28	29	30	31	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49		

AC10P




Shock position: backrest parallel to the floor, seat/ legs inclined depending on the choice of the AC35 series..

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	6	7	9	10	12	14	15	16	17	20	22	27	33	34	35	37	38	40	41	43	45
47	49																						


AC11



Roll holder positioned on the back of the armchair, in stainless steel for rolls with a diameter of 140/600 mm.

COMPATIBLE FOR THE FOLLOWING CHAIRS

- | | | |
|---|---|---|
| 6 | 7 | 8 |
|---|---|---|



Roll holder positioned on the back of the armchair, in stainless steel for rolls with a diameter of 140/600 mm.
With standard armrests.

COMPATIBLE FOR THE FOLLOWING CHAIRS

- | | | |
|---|---|---|
| 1 | 3 | 4 |
|---|---|---|



Roll holder positioned on the back of the armchair, in stainless steel for rolls with a diameter of 140/600 mm.
With multifunction armrests.

COMPATIBLE FOR THE FOLLOWING CHAIRS

- | | | | |
|---|---|---|---|
| 1 | 2 | 3 | 4 |
|---|---|---|---|



Roll holder positioned on the back of the armchair, in stainless steel for rolls with a diameter of 140/600 mm.

COMPATIBLE FOR THE FOLLOWING CHAIRS

- | | | |
|---|----|----|
| 9 | 10 | 11 |
|---|----|----|

AC11/1



Roll holder positioned on the back of the armchair, in stainless steel for rolls with a diameter of 140/600 mm.


After-sales supply.

List price does not include assembly.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	3	4	6	7	8	9	10	11	20	21	22	23	24	27	33	35	36	37	38	39	40	41	42
45	47	49																					

AC11/C



Roll holder positioned on the back of the armchair, in stainless steel for rolls with a diameter of 140/600 mm.

COMPATIBLE FOR THE FOLLOWING CHAIRS

12	13	14	15	16	17	31	43	44	48
----	----	----	----	----	----	----	----	----	----



AC12




Independent movements of the backrest and leg sections. Additional engine to the standard equipment.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	19	20	21	22	23	24	25	26
27	28	29	30	31	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	49			

AC13




Motorized telescopic footrest with HPL shelf and support made of steel, adjustment via pushbutton panel.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	20	21	22	23	27	29	33	34	35	36	37	38	39	40
41	42	45	47	49																			

AC13/1




Motorized telescopic footrest equipped with further manual extension with HPL shelf and support made of steel, adjustment via pushbutton panel

COMPATIBLE FOR THE FOLLOWING CHAIRS

[illegible]



AC13/C1



Motorized footrest made in painted steel. The surface is reversible and therefore adaptable for patients of all heights. Adjustment via pushbutton. Easily removable.
Support surface with padding that can be created at different heights according to needs. The padding is covered with wear-resistant imitation leather material.

COMPATIBLE FOR THE FOLLOWING CHAIRS

12	13	14	43	44
----	----	----	----	----




Support surface with padding that can be created at different heights according to needs. The padding is covered with wear-resistant imitation leather material.

12 13 14 43 44

[illegible]

AC14P



Basket with stainless steel lid with revolving and height-adjustable support, positioned on the accessory rod. Dimensions 250x120x60h mm

On request it is available with different dimensions

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	19	20	21	22	23	24	27	29
31	33	34	35	36	37	38	39	40	41	42	43	44	45	47	48	49							

AC14P/1



Basket with removable and washable **stainless steel lid**, positioned on the armrest. Dimensions 250x120x60h mm

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	19	20	21	22	23	24	26	27
29	31	33	34	35	36	37	38	39	40	41	42	43	44	45	47	48	49						

AC15/PLUS

Serving table with double shelf made of ABS. The frame, on which the swivel wheels that allow movement are installed, is made of painted steel.

The shelf can be tilted into a lectern position and height-adjustable using pneumatic mechanisms. The load capacity is 10 kg.

Dimensions of the first shelf 420x200 mm.
Dimensions of the second shelf 420x550 mm.
Minimum height 830 mm.
Maximum height 1100 mm.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	4	6	8	9	11	12	13	14	15	16	17	19	20	22	24	25	26	27	28	29	30	31	33
36	37	39	40	42	43	44	45	46	47	48	49												

AC15/TOP



Serving table with shelf made of ABS.
The frame, on which the swivel wheels that allow movement are installed, is made of painted steel.

The shelf can be tilted into a lectern position and height-adjustable using pneumatic mechanisms. The load capacity is 10 kg.

Shelf dimensions 420x550 mm.
Minimum height 830 mm.
Maximum height 1100 mm.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	4	6	8	9	11	12	13	14	15	16	17	19	20	22	24	25	26	27	28	29	30	31	33
36	37	39	40	42	43	44	45	46	47	48	49												

AC17P/A



N. 2 batteries with battery charger (4.5 Ah).

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	15	20	21	22	23	24	26	27	29	30	33	34	35	36
37	38	39	40	41	42	45	47	49															

AC19



Armrest mechanism made of stainless materials.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	20	21	22	23	24	33	34	35	36	37	38	39	40	41
42																							

AC21/HP2



Elliptical side walls made of HPL.

Not compatible
with external armrests (AC30H1).

COMPATIBLE FOR THE FOLLOWING CHAIRS

29

AC25P




Support for button panel housing integrated under the armrest of the chair..

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	19	20	21	22	23	24	26	27
29	30	31	33	34	35	36	37	38	39	40	41	42	43	44	45	47	49						

AC26



Weight loss alarms.

COMPATIBLE FOR THE FOLLOWING CHAIRS

33	34	35	36	37	38	39	40	41	42	43	44	45	46	47
----	----	----	----	----	----	----	----	----	----	----	----	----	----	----

AC27



Additional headrest cushion adjustable in height, with padding made of expanded resin and antibacterial imitation leather covering.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	19	20	21	22	23	24	26	27
31	33	34	35	36	37	38	39	40	41	42	43	44	45	47	48	49							

AC27/1



Additional headrest cushion adjustable in height, with padding made of expanded resin and antibacterial imitation leather covering.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	19	20	21	22	23	24	26	27
31	33	34	35	36	37	38	39	40	41	42	43	44	45	47	48	49							

AC27/2



Additional headrest cushion adjustable in height, with padding made of expanded resin and antibacterial imitation leather covering.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	19	20	21	22	23	24	26	27
31	33	34	35	36	37	38	39	40	41	42	43	44	45	47	48	49							

AC27/3



Additional headrest cushion radjustable in height, with padding made of expanded resin and antibacterial imitation leather covering.

MEDIUM DENSITY

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	19	20	21	22	23	24	26	27
31	33	34	35	36	37	38	39	40	41	42	43	44	45	47	48	49							

AC27/4




Additional headrest cushion adjustable in height, with padding made of expanded resin and antibacterial imitation leather covering.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	19	20	21	22	23	24	26	27
31	33	34	35	36	37	38	39	40	41	42	43	44	45	47	48	49							

AC27/5



Additional headrest cushion adjustable in height, with padding made of expanded resin and antibacterial imitation leather covering.

COMPATIBLE FOR THE FOLLOWING CHAIRS

29

AC27/6




Additional headrest cushion adjustable in height, with padding made of expanded resin and antibacterial imitation leather covering.

LOW DENSITY

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	19	20	21	22	23	24	26	27
31	33	34	35	36	37	38	39	40	41	42	43	44	45	47	48	49							

AC27/7



Additional headrest cushion adjustable in height, with padding made of expanded resin and antibacterial imitation leather covering.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	19	20	21	22	23	24	26	27
31	33	34	35	36	37	38	39	40	41	42	43	44	45	47	48	49							

AC27/8



Additional headrest cushion adjustable in height, with padding made of expanded resin and antibacterial imitation leather covering.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	19	20	21	22	23	24	26	27
31	33	34	35	36	37	38	39	40	41	42	43	44	45	47	48	49							

AC27/9



Additional headrest cushion adjustable in height, with padding made of expanded resin and antibacterial imitation leather covering.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	19	20	21	22	23	24	26	27
31	33	34	35	36	37	38	39	40	41	42	43	44	45	47	48	49							

AC27/10



Additional headrest cushion adjustable in height, with padding made of expanded resin and antibacterial imitation leather covering.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	19	20	21	22	23	24	26	27
31	33	34	35	36	37	38	39	40	41	42	43	44	45	47	48	49							

AC28/200-248



Variable height adjustment via motorized telescopic column with rectangular section with a capacity of 200 kg and an excursion of 248 mm

COMPATIBLE FOR THE FOLLOWING CHAIRS

24

AC28/200-400



Variable height adjustment via motorized telescopic column with rectangular section with a capacity of 200 kg and an excursion of 400 mm.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	6	9	12	14	15	16	17	19	27	31	33	37	40	43	45	47	49
---	---	---	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----

AC28/300



Variable height adjustment via motorized telescopic column with rectangular section with a capacity of 300 kg and an excursion of 250 mm. For bariatric chairs.

COMPATIBLE FOR THE FOLLOWING CHAIRS

4	8	11	13	29	36	39	42	44
---	---	----	----	----	----	----	----	----

AC29




Armrest band made of imitation leather, adjustable via Velcro.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	19	20	21	22	23	24	26	27
28	29	30	31	33	34	35	36	37	38	39	40	41	42	43	44	45	48						

AC30P




Multifunction armrest in stainless steel and aluminum adjustable in: height, adjustable, tilting, folding and removable.

COMPATIBLE FOR THE FOLLOWING CHAIRS

[illegible]

AC30/C1



Ergonomic multifunction armrest, adjustable via knob.

Movements allowed:

- integral with the backrest (foto 1)
- swivel (foto 2)
- adjustable in height (foto 3)
- removable (foto 4)
- foldable (foto 5)
- tiltable “up/down” (foto 6).

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	19	20	21	22	23	24	27	31
33	34	35	36	37	38	39	40	41	42	43	44	45	48										


Implementation of the AC30/C1 armrest with the **function tion slide**.

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	19	20	21	22	23	24	27	31
33	34	35	36	37	38	39	40	41	42	43	44	45	48										

Implementation of the AC30/C1 armrest with the **internal rotation function.**

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	19	20	21	22	23	24	27	31
33	34	35	36	37	38	39	40	41	42	43	44	45	48										

AC30/C5




Slim multifunction armrest in steel adjustable via knob with slide function.
Armrest dimensions: length 515 mm, width 105 mm..

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	19	20	21	22	23	24	27	31
33	34	35	36	37	38	39	40	41	42	43	44	45	48										

AC30/C6



Adjustable multifunctional armrest with gas spring.

COMPATIBLE FOR THE FOLLOWING CHAIRS

27	45	47	49
----	----	----	----

AC30/H1



Multifunctional armrest adjustable in height, swivel and removable if necessary.

COMPATIBLE FOR THE FOLLOWING CHAIRS

29

AC31



Joint cushion rmade with expanded resin and antibacterial imitation leather covering. Slim multifunction armrest in steel adjustable via knob with slide function. Thanks to its ergonomics and its cylindrical shape, it is functional to optimize comfort during long stays..

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	20	21	22	23	27	33	34	35	36	37	38	39	45
---	---	---	---	---	---	---	----	----	----	----	----	----	----	----	----	----	----	----	----

AC32P



Courtesy **LED light** equipped with “ON/OFF” easily accessible to the patient

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	20	21	22	23	27	29	33	34
35	36	37	38	39	40	41	42	43	44	45	47	49											

AC33



Manually **height-adjustable headrest** with 4 intermediate positions. Mechanism made of stainless material and easily removable if necessary.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	19	20	21	22	23	27	31	33	34	35	36	37	38	39	45	48
---	---	---	---	---	---	---	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----

AC33/1



Reclining headrest 0/15°.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	4	6	8	33	36	37	39	47	49
---	---	---	---	----	----	----	----	----	----

AC33/2



Headrest reclining 0/15° and adjustable in size manually.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	4	6	8	33	36	37	39	47	49
---	---	---	---	----	----	----	----	----	----

AC34P



Key for inhibiting push button controls. This device, supplied to healthcare personnel, can be used to protect patient safety.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	19	20	21	22	23	24	25	26
27	28	29	30	31	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	49			

AC35



Fixed seat parallel to the floor.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	33	34	35	36
---	---	---	---	----	----	----	----

AC35/1



Motorized seat with adjustable inclination from 0° to -15°, operated via pushbutton. This function is useful for passive postural variation of the patient for the prevention of bedsores.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	19	20	22	24	26	27	29	31
33	34	35	36	37	38	39	40	41	42	43	44	45	47	49									

AC35/2



Motorized seat with adjustable inclination from +3° to -12°, operated via pushbutton. This function is useful for passive postural variation of the patient for the prevention of bedsores..

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	6	7	9	10	20	22	24	27	33	34	35	37	38	40	41	45	47	49
---	---	---	---	---	---	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----

AC35/3



Motorized seat with adjustable inclination from -10° to -25°, operated via pushbutton. This function is useful for passive postural variation of the patient for the prevention of bedsores.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	6	7	9	10	20	22	24	27	33	34	35	37	38	40	41	45	47	49
---	---	---	---	---	---	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----

AC35/4



Motorized seat with adjustable inclination from -12° to +12°, operated via pushbutton.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	33	34	35	47	49
---	---	---	----	----	----	----	----

AC35/5



Motorized seat with adjustable inclination from 0° to +10°, operated via buttons.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	33	34	35	47	49
---	---	---	----	----	----	----	----


AC36



Pair of removable transparent PVC covers to protect the armrests, useful for preserving the integrity of the upholstery over time.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	19	20	21	22	23	24	27	28
29	30	31	33	34	35	36	37	38	39	40	41	42	43	44	45	47	48	49					



Pair of fixed transparent PVC covers to protect the armrests, useful for preserving the integrity of the upholstery over time.


1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	19	20	21	22	23	24	27	28
29	30	31	33	34	35	36	37	38	39	40	41	42	43	44	45	47	48	49					

A close-up photograph of a teal-colored car seat and headrest. The seat is on the left, and the headrest is on the right, both featuring a smooth, glossy finish. The background is a plain, light gray.

Transparent PVC covers for footrests, useful for preserving the integrity of the upholstery over time.

1	2	3	4	6	7	8	9	10	11	12	13	14	20	21	22	23	24	27	28	29	30	33	34
35	36	37	38	39	40	41	42	43	44	45	47	49											


AC36/2



Transparent PVC cover for the leg section
integral with the cushion, useful for preserving the integrity of the upholstery over time.
On request it can also be made removable.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	20	21	22	23	28
30	32	33	34	35	36	37	38	39			



Transparent PVC cover for the leg section
integral with the cushion, useful for preserving the integrity of the upholstery over time.
On request it can also be made removable.

COMPATIBLE FOR THE FOLLOWING CHAIRS

12	13	14	15	16	17	18	19	31	43	44	48
----	----	----	----	----	----	----	----	----	----	----	----



Transparent PVC cover for the leg section
integral with the cushion, useful for preserving the integrity of the upholstery over time.
On request it can also be made removable.

COMPATIBLE FOR THE FOLLOWING CHAIRS

[illegible]

AC37



Extra large armrest 200 mm..

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	20	21	22	23	24	27	33	34
35	36	37	38	39	40	41	42	43	44	45													

AC37/1



Bracciolo extra long 625x140 mm.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	20	21	22	23	27	33	34	35	36	37	38	39	40	41
42	45																						

AC37/3



Elongated armrest made with expanded resin on the back.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	20	22	27	33	34	35	36	37	38	39	40	41	42	45
---	---	---	---	---	---	---	---	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----

AC37/4



Seamless armrest supplement.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	20	21	22	23	27	33	34	35
36	37	38	39	40	41	42	43	44	45														

AC37/A



Extra large armrest 360x120 mm.

COMPATIBLE FOR THE FOLLOWING CHAIRS

32

AC38



Additional cushion for lowered seat:

- Padding with standard height (photo 1).
- Seat in the lowered version (photo 2).

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	20	21	22	23	27	33	34	35	36	37	38	39	45
---	---	---	---	---	---	---	----	----	----	----	----	----	----	----	----	----	----	----	----

AC39



Leg section at 90° to the floor.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	27	33	34	35	36	37	38	39	45
---	---	---	---	---	---	---	----	----	----	----	----	----	----	----	----

A close-up photograph of a red leather examination table. The table is positioned at an angle, showing the backrest and the main table surface. A manual air pump with a black hose is attached to the side of the table, likely for inflating an air cushion. A metal control lever is visible near the pump. The table is set against a plain white background.


Air-adjustable internal lumbar inflatable cushion. Functional to optimize posture during the phases of therapy with preventive antidecubitus effect.

1	2	3	4	6	7	8	9	10	11	12	13	14	20	21	22	23	27	33	34	35	36	37	38
39	40	41	42	43	44	45	47	49															

A purple inflatable cushion, likely for a car seat, with a black pump and hose attached to the top. The cushion is rectangular and appears to be made of a smooth, possibly vinyl or plastic material. The pump is a small black device with a hose that loops around the top of the cushion.

Air-adjustable external lumbar inflatable cushion. Functional to optimize posture during the phases of therapy with preventive antidecubitus effect.

1	2	3	4	6	7	8	9	10	11	12	13	14	20	21	22	23	24	26	27	28	29	30	33
34	35	36	37	38	39	40	41	42	43	44	45	47	49										



Additional external lumbar cushion. Functional to optimize posture during the phases of therapy with a preventive anti-decubitus effect.

1	2	3	4	6	7	8	9	10	11	12	13	14	20	21	22	23	24	26	27	28	29	30	33
34	35	36	37	38	39	40	41	42	43	44	45	47	49										

AC41



Thermal printer for armchairs with weighing system, useful for printing the weight value in real time.

It is powered by rechargeable batteries and the supplied power supply, and can also be fixed to the chair with a dedicated support..

COMPATIBLE FOR THE FOLLOWING CHAIRS

33	34	35	36	37	38	39	40	41	42	43	44	45	46	47
----	----	----	----	----	----	----	----	----	----	----	----	----	----	----

AC43P




Sick lift pole made of reinforced steel tubing and ergonomic triangular handle, adjustable in height via a special mechanism.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	4	6	8	9	11	20	22	27	29	33	36	37	39	40	42	45	47	49
---	---	---	---	---	----	----	----	----	----	----	----	----	----	----	----	----	----	----

AC48




TV 24" installed on an armchair using a dedicated swivel support.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	12	13	14	20	21	22	23	27	33	34	35	36	37	38
39	40	41	42	43	44	45	47	49															

AC48/2



TV 24" installed on armchair via dedicated swivel mount.

COMPATIBILE PER LE SEGUENTI POLTRONE

1	2	3	4	6	7	8	9	10	11	12	13	14	20	21	22	23	27	33	34	35	36	37	38
39	40	41	42	43	44	45	47	49															

AC48/PP



Preparation for TV support..

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	12	13	14	20	21	22	23	27	33	34	35	36	37	38
39	40	41	42	43	44	45	47	49															

AC48/SP



TV stand.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	12	13	14	20	21	22	23	27	33	34	35	36	37	38
39	40	41	42	43	44	45	47	49															

Pair of side panels made of HPL, functional to preserving the integrity of the base as well as performing its aesthetic function.

1	2	3	6	7	9	10	16	17	27	33	34	35	37	38	40	41	45
---	---	---	---	---	---	----	----	----	----	----	----	----	----	----	----	----	----



Cable winding support

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	20	21	22	23	24	27	29	30
31	33	34	35	36	37	38	39	40	41	42	43	44	45	47	49								

Backrest made of memory foam.

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	20	21	22	23	27	33	34	35
36	37	38	39	40	41	42	43	44	45	47	49												

AC52



Seat made of memory foam

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	20	21	22	23	24	27	33	34
35	36	37	38	39	40	41	42	43	44	45	47	49											

AC55P



Built-in servant’s table removable if necessary. Dimension 420x550 mm.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	4	6	8	9	11	12	13	14	20	22	27	33	36	37	39	40	42	43	44	45
---	---	---	---	---	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----

AC56




Handle made of stainless steel to facilitate transport activities.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	4	6	8	9	11	20	22	24	27	28	29	30	33	36	37	39	40	42	45	47	49
---	---	---	---	---	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----

AC56/1




Pvc transport handles.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	20	21	22	23	27	28	30	33
34	35	36	37	38	39	40	41	42	43	44	45	47	49										

AC56/C



Handle made of stainless steel to facilitate transport activities.

COMPATIBLE FOR THE FOLLOWING CHAIRS

12	13	14	15	16	17	19	43	44	48
----	----	----	----	----	----	----	----	----	----

AC56/C1



Handle made of aluminum to facilitate transport activities.

COMPATIBLE FOR THE FOLLOWING CHAIRS

12	13	14	15	16	17	31	43	44
----	----	----	----	----	----	----	----	----

AC57



Print for AC49 side panels to protect the base.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	6	7	9	10	27	33	34	35	37	38	40	41	45
---	---	---	---	---	---	----	----	----	----	----	----	----	----	----	----

AC60



Adjustable belts to protect the patient's safety

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	20	21	22	23	27	29	30	33
34	35	36	37	38	39	40	41	42	43	44	45	47	49										

AC63




Leg section at 90° to the seat

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	6	7	16	17	19	20	22	27	28	31	33	34	35	37	38	45	48
---	---	---	---	---	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----

AC63/1



Leg section at 90° to the seat for bariatric chairs.

COMPATIBLE FOR THE FOLLOWING CHAIRS

4	8	36	39
---	---	----	----

Automatic emergency button to obtain Trendelenburg or shock positions.

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	20	22	24	27	29	33	34	35
36	37	38	39	40	41	42	43	44	45	47	49												

A close-up photograph of a black PlayStation 4 console. A green circular sticker with a stylized 'i' logo is affixed to the front. The console's disc cover is partially visible, showing a green center. A small green light bar is visible below the disc cover.

Button to obtain automatic input position.

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	19	20	22	27	31	33	34	35
36	37	38	39	40	41	42	43	44	45	47	49												

Button to obtain the clinostati position.

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	20	22	27	29	33	34	35	36
37	38	39	40	41	42	43	44	45	47	49													

A 3D rendering of a red, glossy, cylindrical object, possibly a button or a small container, sitting on a white circular base. The object has a smooth, reflective surface and is centered on a plain white background.

Emergency stop button.

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	19	20	22	26	27	29	31	33
34	35	36	37	38	39	40	41	42	43	44	45	47	49										



Adjustable cup holder located under the armrest.

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	19	20	21	22	23	24	27	31
33	34	35	36	37	38	39	40	41	42	43	44	45	48										

AC75/1



Adjustable courtesy table placed under the armrest with slots for glasses or smartphones.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	4	6	8	11	12	13	19	20	22	24	27	31	33	36	37	39	42	43	44	45	48
---	---	---	---	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----

AC75/2





Adjustable cup holder placed under the armrest, extra-long version..

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	4	6	8	11	12	13	19	20	22	24	27	31	33	36	37	39	42	43	44	45	48
---	---	---	---	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----

Multifunctional adjustable support with hooks located under the armrest, useful for holding small bags.



1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	20	21	22	23	24	27	31	33
34	35	36	37	38	39	40	41	42	43	44	45	48											

Additional cushion for back section.



1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	20	21	22	23	24	27	33	34
35	36	37	38	39	40	41	42	43	44	45	47	49											

AC79



Built-in thermo-cushion for backrest or seat.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	20	21	22	23	27	33	34	35
36	37	38	39	40	41	42	43	44	45	47	49												

AC81P



Articulated and removable servant's table, in stainless steel and aluminum with storage edge. Made of non-porous material without bacterial retention.

Top dimensions equal to 440x340 mm.



COMPATIBLE FOR THE FOLLOWING CHAIRS

1	4	6	8	9	11	12	14	20	22	27	33	36	37	39	40	42	43	44	45	47	49
---	---	---	---	---	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----

AC81/SP




Arrangement for serving table accommodation.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	4	6	8	9	11	12	14	20	22	27	33	36	37	39	40	42	43	44	45	47	49
---	---	---	---	---	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----

AC82



Leg section with motorized telescopic extension, useful in the case of configuration with 90° leg section to give greater comfort to the patient.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	6	27	33	37	45
---	---	----	----	----	----

AC82/1




Leg section with motorized telescopic extension for bariatric chairs.

COMPATIBLE FOR THE FOLLOWING CHAIRS

4	8	36	39
---	---	----	----

AC83




Cover made of ABS for base. Functional to protect and preserve the sensitive parts of the device.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	6	9	12	13	14	15	26	27	33	37	40	43	45	47	49
---	---	---	----	----	----	----	----	----	----	----	----	----	----	----	----

AC83/C




Cover made of ABS for base. Functional to protect and preserve the sensitive parts of the device.

COMPATIBLE FOR THE FOLLOWING CHAIRS

16	17	19	31	48
----	----	----	----	----

AC83/H



Cover made of ABS for base. Functional to protect and preserve the sensitive parts of the device.

COMPATIBLE FOR THE FOLLOWING CHAIRS

4	8	11	13	29	44
---	---	----	----	----	----

AC84



Cover made of ABS for seat. Functional to protect and preserve the sensitive parts of the device.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	6	9	20	22	27	33	37	40	45	47	49
---	---	---	----	----	----	----	----	----	----	----	----

AC84/C



Cover made of ABS for seat section. Functional to protect and preserve the sensitive parts of the device.

COMPATIBLE FOR THE FOLLOWING CHAIRS

12	14	15	16	17	19	31	43	48
----	----	----	----	----	----	----	----	----

AC84/H



Cover made of ABS for seat section. Functional to protect and preserve the sensitive parts of the device.

COMPATIBLE FOR THE FOLLOWING CHAIRS

29

AC85/C



Cover made of ABS for leg section. Functional to protect and preserve the sensitive parts of the device.

COMPATIBLE FOR THE FOLLOWING CHAIRS

- | | | | |
|----|----|----|----|
| 12 | 14 | 15 | 43 |
|----|----|----|----|

AC86/C



Cover made of ABS for backrest section. Functional to protect and preserve the sensitive parts of the device.

COMPATIBLE FOR THE FOLLOWING CHAIRS

- | | | | | | | | |
|----|----|----|----|----|----|----|----|
| 12 | 14 | 15 | 16 | 17 | 19 | 31 | 43 |
|----|----|----|----|----|----|----|----|

AC86/H



Cover made of ABS for backrest section. Functional to protect and preserve the sensitive parts of the device.

COMPATIBLE FOR THE FOLLOWING CHAIRS

- | |
|----|
| 29 |
|----|

AC87



Bi-elastic silicone coating, extremely resistant and stain-resistant. To choose the color, consult the appropriate color card.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	19	20	21	22	23	24	25	26
27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	

AC87/1



Leatherette covering with stain-resistant properties. To choose the color, consult the appropriate color card.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	19	20	21	22	23	24	25	26
27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	

AC88



Armrest cover made of material: waterproof, washable, breathable, bi-elastic and bacteriostatic.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	19	20	21	22	23	24	27	28
29	30	31	33	34	35	36	37	38	39	40	41	42	43	44	45	48							

AC89



Armchair cover made of material: waterproof, washable, breathable, bi-elastic and bacteriostatic.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	19	20	21	22	23	24	26	27
28	29	30	31	33	34	35	36	37	38	39	40	41	42	43	44	45	47	48	49				

AC90/H



Armchair cover in TNT material.

COMPATIBLE FOR THE FOLLOWING CHAIRS

29

AC90/H1



Cotton **armchair cover**.

COMPATIBLE FOR THE FOLLOWING CHAIRS

29

AC94P




Wireless call kit installed on the armrest of the chair.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	20	21	22	23	24	27	28	29
30	33	34	35	36	37	38	39	40	41	42	43	44	45										

AC95P



Cylinder holder installed on the base of the chair.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	4	6	8	9	11	12	13	14	15	16	17	20	22	26	27	33	36	37	39	40	42	43	44
45	47	49																					

AC96P



Wireless call button installed on the armrest of the chair.

COMPATIBLE FOR THE FOLLOWING CHAIRS

26

AC99



90° backrest section.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	20	21	22	23	27	33	34	35
36	37	38	39	40	41	42	43	44	45	47	49												

AC109



Adjustable side rails installed on DIN bar.

Incompatible
with the Stephen H Comfort 3.0 and Stephen Scale H Comfort 3.0.

COMPATIBLE FOR THE FOLLOWING CHAIRS

6	8	12	13	26	28	37	39	43	44
---	---	----	----	----	----	----	----	----	----

AC109/1



Adjustable side rails installed on DIN bar.

Incompatible
with the Stephen H Comfort 3.0 and Stephen Scale H Comfort 3.0.

COMPATIBLE FOR THE FOLLOWING CHAIRS

6	8	12	22	26	28	37	39	43
---	---	----	----	----	----	----	----	----

AC110



Leg rests installed on a DIN bar adjustable using a special mechanism.

COMPATIBLE FOR THE FOLLOWING CHAIRS

6	8	22	37	39
---	---	----	----	----

AC115




Activation of the controls for moving the sections on the Fly Plus display:
backrest, legs, variable height and memory positions:
entrance position and emergency position.

COMPATIBLE FOR THE FOLLOWING CHAIRS

33	34	35	36	37	38	39	40	41	42	43	44	45	46	47
----	----	----	----	----	----	----	----	----	----	----	----	----	----	----

AC116



Flexible document holder positioned on the back of the document chair. Size 260x350 mm.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	19	20	21	22	23	24	25	26
27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	

AC117



Multiple transport handle for Athe-na 2.0,
useful for transport both in a postural position
and in a stretcher position. These handles
facilitate maneuverability of the device even with
the patient on board.

COMPATIBLE FOR THE FOLLOWING CHAIRS

26

AC118



5^a directional wheel with pedal activation to facilitate movement for long journeys from place to place.

COMPATIBLE FOR THE FOLLOWING CHAIRS

26

AC119



Pelvic girdle.

COMPATIBLE FOR THE FOLLOWING CHAIRS

12 13 14 15 16 17 26 43 44

AC122



Independent pneumatic movement of the backrest and legrest sections.

COMPATIBLE FOR THE FOLLOWING CHAIRS

32

Compass sides with servo-assisted folding mechanism made of polymer and aluminium.

26

Seat width extendable by 20 cm on DIN rail.

1	4	6	8	9	11	12	13	14	15	16	20	22	26	27	29	33	36	37	39	40	42	43	44
---	---	---	---	---	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----

45


Bilateral steel DIN bar for housing additional accessories measuring 25x8 mm.

1	4	6	8	9	11	12	13	14	15	16	20	24	22	26	27	29	33	36	37	39	40	42	43
---	---	---	---	---	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----

44 45 47 49

AC135

Support feet with non-slip base.

COMPATIBLE FOR THE FOLLOWING CHAIRS**32****AC136**

N. 2 rear directional wheels diameter 125 mm equipped with anti-slip front tips.

COMPATIBLE FOR THE FOLLOWING CHAIRS**32****AC144**

Multifunction armrest with radial joint installed on DIN bar.

COMPATIBLE FOR THE FOLLOWING CHAIRS**26**

AC145



Support with radial joint on DIN bar and IV pole.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	4	6	8	9	11	12	13	14	15	16	20	22	26	27	29	33	36	37	39	40	42	43	44
45	47	49																					

AC146



USB port to allow charging/powering of electronic devices with safety cover.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	20	21	22	23	24	25	26	27
28	29	30	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	49				

AC150



Leg section equipped with footrest with folding HPL top.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	4	6	8	20	22	27	33	36	37	39	45
---	---	---	---	----	----	----	----	----	----	----	----

AC153



Footrest with swivel mechanism and automatic braking system.

COMPATIBLE FOR THE FOLLOWING CHAIRS

- 16
- 17

AC157




Interchangeable articulated headrest

COMPATIBLE FOR THE FOLLOWING CHAIRS

- 15
- 17

AC159



Column stand for battery chargers.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	19	20	21	22	23	24	25	26
27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	

AC164



90° rotatable display for scale chairs.

COMPATIBLE FOR THE FOLLOWING CHAIRS

- | | | | | | | | | | | | | | | |
|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|
| 33 | 34 | 35 | 36 | 37 | 38 | 39 | 40 | 41 | 42 | 43 | 44 | 45 | 46 | 47 |
|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|

AC167



GBTouch multifunction control touch screen display.

COMPATIBLE FOR THE FOLLOWING CHAIRS

- | | | | | | | | | | | | | | | |
|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|
| 33 | 34 | 35 | 36 | 37 | 38 | 39 | 40 | 41 | 42 | 43 | 44 | 45 | 46 | 47 |
|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|

AC167/1




WI-FI/Bluetooth combo module for GBTouch.

COMPATIBLE FOR THE FOLLOWING CHAIRS

- | | | | | | | | | | | | | | | |
|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|
| 33 | 34 | 35 | 36 | 37 | 38 | 39 | 40 | 41 | 42 | 43 | 44 | 45 | 46 | 47 |
|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|

AC167/2




GBTouch **USB stick for storage**

COMPATIBLE FOR THE FOLLOWING CHAIRS

33	34	35	36	37	38	39	40	41	42	43	44	45	46	47
----	----	----	----	----	----	----	----	----	----	----	----	----	----	----

AC167/3




10/100 ethernet port for GBTouch.

COMPATIBLE FOR THE FOLLOWING CHAIRS

33	34	35	36	37	38	39	40	41	42	43	44	45	46	47
----	----	----	----	----	----	----	----	----	----	----	----	----	----	----

AC167/4




USB A port for GBTouch. Weight data transmissions.
Not compatible with the AC167/2 USB stick.

COMPATIBLE FOR THE FOLLOWING CHAIRS

33	34	35	36	37	38	39	40	41	42	43	44	45	46	47
----	----	----	----	----	----	----	----	----	----	----	----	----	----	----

AC167/5




USB B port for GBTouch. External printer connection.

COMPATIBLE FOR THE FOLLOWING CHAIRS

- | | | | | | | | | | | | | | | |
|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|
| 33 | 34 | 35 | 36 | 37 | 38 | 39 | 40 | 41 | 42 | 43 | 44 | 45 | 46 | 47 |
|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|

AC167/6




RS-232 port for GBTouch.

COMPATIBLE FOR THE FOLLOWING CHAIRS

- | | | | | | | | | | | | | | | |
|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|
| 33 | 34 | 35 | 36 | 37 | 38 | 39 | 40 | 41 | 42 | 43 | 44 | 45 | 46 | 47 |
|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|

AC167/7



RS-485 port for GBTouch.

COMPATIBLE FOR THE FOLLOWING CHAIRS

- | | | | | | | | | | | | | | | |
|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|
| 33 | 34 | 35 | 36 | 37 | 38 | 39 | 40 | 41 | 42 | 43 | 44 | 45 | 46 | 47 |
|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|

AC174



Legrest with leg section 90°/180° with respect to the floor and with retractable folding footrest.

COMPATIBLE FOR THE FOLLOWING CHAIRS

- | | | | |
|---|---|----|----|
| 1 | 6 | 33 | 37 |
|---|---|----|----|

AC186



Shaped back section for armrest AC30/C6.

COMPATIBLE FOR THE FOLLOWING CHAIRS

- | | |
|----|----|
| 47 | 49 |
|----|----|

Once the device has been configured, consult the Gardhen Bilance sales offices in order to avoid any incompatibilities between the chosen accessories.

The images reproduced are purely indicative and the manufacturer has the right to make changes to the projects without notice.



Gardhen Bilance S.r.l.

Via G. Luraghi c/o Consorzio Il Sole

lotto S - 80038 Pomigliano d'Arco (NA) - Italy

Tel. +39 081 8692160 (through selection 5 lines)

+39 3913087453 (if the fixed network is not reachable)

E-mail: info@gardhenbilance.it

PEC: gardhenbilance@pec.it