


Positioner Wedge




**SINGLE USE
POSITIONER**

Frequent use


Surgery

Emergency room

Intensive care unit



PRESSURE REDUCTION AREA SCHEME



An important phase of any surgery is the accommodation of the patient on the operating table. In fact correct positioning reduces the risk of ulcer formation, pressure sores and nerve injuries during the operation.

Furthermore, the patient who goes through local or total anesthesia gradually loses consciousness, motor skills or at least sensitivity to touch and pain. He is therefore unable to perceive and react to potentially harmful external stimuli related to the position in which he is placed on the operating table or caused involuntarily and unknowingly by the health and social care staff: friction, rubbing, excessive and prolonged pressure, shocks, burns.

Therefore it is advisable to choose appropriate positioning devices that ensure correct pressure distribution and maintain the patient in the proper position throughout the intervention.

Reliable, efficient and not much expensive, the positioners are specially designed tools to prevent the onset of injuries and to relieve and stabilize any injuries already present in the patient undergoing surgery.

The wedge positioner can be used for positioning and support of the back and lower limbs in different types of operations.

The colour of the foam used is subject to availability

Benefits



- Effectiveness to prevent abrasions;
- Effectiveness in preventing friction burns;
- Effectiveness to prevent decubitus wounds;
- Ensures air circulation;
- Ergonomic design.

Features



- Polyurethan foam;
- 100% latex free;
- Non-sterile;
- Single-use.

Dimensions



| ID Product | Dimensions |
|--------------|----------------------------|
| GB.0006PUM01 | L. 20 x P. 46 x H. 20 cm |
| GB.0006PUM02 | L. 31 x P. 56 x H. 23 cm |
| GB.0006PUM03 | L. 60 x P. 58 x H. 22,5 cm |
| GB.0006PUM04 | L. 45 x P. 60 x H. 45 cm |
| GB.0006PUM05 | L. 30 x P. 55 x H. 22 cm |

Technical sheet



| | |
|--------------|-----------------------|
| GMDN Code | 41058 |
| Basic UDI-DI | 805771740POSPU |
| Intended use | Upholstered cushion |
| Manufacturer | GARDHEN BILANCE SRL |
| Class | Class I |

For certifications please refer to "Certifications and Quality system" section

Pictures contained in this brochure is for illustrative purpose only. The manufacturer reserves the right to modify the projects without notice. For accessories identification please refer to the accessory catalogue or commercial offer. When data are different please refer to commercial offer.

Dep. EN_255 rev. 06/05/2022

