





An important phase of any surgery is the accommodation of the patient on the operating table. In fact correct positioning reduces the risk of ulcer formation, pressure sores and nerve injuries during the operation.

Furthermore, the patient who goes through local or total anesthesia gradually loses consciousness, motor skills or at least sensitivity to touch and pain. He is therefore unable to perceive and react to potentially harmful external stimuli related to the position in which he is placed on the operating table or caused involuntarily and unknowingly by the health and social care staff: friction, rubbing, excessive and prolonged pressure, shocks, burns.

Therefore it is advisable to choose appropriate positioning devices that ensure correct pressure distribution and maintain the patient in the proper position troughout the intervention.

Reliable, efficient and not much expensive, the positioners are specially designed tools to prevent the onset of injuries and to relieve and stabilize any injuries already present in the patient undergoing surgery.

Gel positioners can be easily cleaned with common used detergents for opthalmic, oral and plastic surgery.

The horseshoe gel positioner ensures the passage of air to the patient in the lateral or prone position.

The color of the gel used for the realization is subject to availability.

Benefits

- Effectiveness to prevent abrasions;
- Effectiveness to prevent friction burns;
- Effectiveness to prevent decubitus wounds;
- Ergonomic design;
- Sanitizable and reusable.

Features

- Polyurethane gel;
- Polyurethane coating;
- Latex free;
- Antibacterial.

Dimensions

ID Product	Basic UDI-DI	Dimensions
GBSK0004	08809240020562	Diam.210 x H 45 mm
GBSK0005	08809240020579	Diam.150 x H 35 mm
GBSK0006	08809240020586	Diam.92 x H 24 mm

Technical sheet

GMDN Code	35229
Intended use	Gelpositioners
Manufacturer	GARDHEN BILANCE SRL
Class	Class I

For certifications please refer to "Certifications and Quality system" section

Pictures contained in this brochure is for illustrative purpose only. The manufacturer reserves the right to modify the projects without notice. For accessories identification please refer to the accessory catalogue or commercial offer. When data are different please refer to commercial offer.

Dep. EN 280 rev. 26/05/2022



2