

Gel Half-cylindrical Positioner



GEL POSITIONERS

Frequent use


Emergency room

Surgery

Intensive care unit



PRESSURE REDUCTION AREA SCHEME



An important phase of any surgery is the accommodation of the patient on the operating table. In fact correct positioning reduces the risk of ulcer formation, pressure sores and nerve injuries during the operation.

Furthermore, the patient who goes through local or total anesthesia gradually loses consciousness, motor skills or at least sensitivity to touch and pain. He is therefore unable to perceive and react to potentially harmful external stimuli related to the position in which he is placed on the operating table or caused involuntarily and unknowingly by the health and social care staff: friction, rubbing, excessive and prolonged pressure, shocks, burns.

Therefore it is advisable to choose appropriate positioning devices that ensure correct pressure distribution and maintain the patient in the proper position throughout the intervention.

Reliable, efficient and not much expensive, the positioners are specially designed tools to prevent the onset of injuries and to relieve and stabilize any injuries already present in the patient undergoing surgery.

Gel positioners can be easily cleaned with common used detergents for ophthalmic, oral and plastic surgery.

The semi-cylindrical gel positioner can be used for positioning and supporting various parts of the body (neck, lower limbs, torso, etc.) in different positions and types of operation.

The colour of the gel in which the positioner is made may change according to availability.

Benefits

- Effectiveness to prevent abrasions;
- Effectiveness to prevent friction burns;
- Effectiveness to prevent decubitus wounds;
- Ergonomic design;
- Sanitizable and reusable.

Features

- Polyurethane gel;
- Polyurethane coating;
- Latex free;
- Antibacterial.

Dimensions

ID Product	Basic UDI-DI	Dimensions
GBSK0034	08809240020821	L.470 x P.100 x H.50 mm
GBSK0035	08809240028803	L.470 x P.150 x H.80 mm
GBSK0036	08809240020883	L.470 x P.150 x H.120 mm
GBSK0037	08809240020890	L.250 x P.150 x H.80 mm
GBSK0038	08809240020906	L.520 x P.100 x H.80 mm

Technical sheet

GMDN Code	35229
Intended use	Gel positioners
Manufacturer	GARDHEN BILANCE SRL
Class	Class I

For certifications please refer to "Certifications and Quality system" section

Pictures contained in this brochure is for illustrative purpose only. The manufacturer reserves the right to modify the projects without notice. For accessories identification please refer to the accessory catalogue or commercial offer. When data are different please refer to commercial offer.

Dep. EN 294 rev. 26/05/2022

