

Securtop *light*

Monobloc countertop

ALL CLINICAL
AREAS

Frequent use

All clinical areas
Psychiatric wards
Mental health facilities



CONTAINMENT & SECURITY

a single structure, maximum strength




www.gardenbalance.it

Securtop^{light}


The **monobloc HPL countertop** is designed to provide durability, hygiene and safety in healthcare settings. Thanks to its structure and direct floor mounting, it is highly stable and resistant to heavy use, minimising maintenance requirements.

The HPL surface is resistant to impacts, wear and tear, and the cleaning agents and disinfectants commonly used in healthcare facilities, making it easy to clean effectively. Its minimalist design, free of superfluous elements, helps to maintain high hygiene standards and facilitates space management. Ideal for hospital wards and mental healthcare areas, this shelf is a reliable, durable solution that combines functionality, safety and comfort.



Standard equipment

- Reinforced frame made of S235JR steel made stainless by pre-treatment and thermosetting epoxy paint;
- Monobloc HPL structure;
- Floor and wall mounting brackets.




Data sheet



CND/EMDN	V0899
GMDN Code	37075
N. progressive/R	3036596/R
UDI-DI di Base	805771740VARPE
Product Code	GB0458.AC221
Intended use	All clinical areas
Manufacturer	GARDHEN BILANCE SRL
Total weight	30 kg
Maximum load evenly distributed	50 kg



Sizes





Data sheet

CND/EMDN	V0899
GMDN Code	37075
N. progressive/R	3036596/R
UDI-DI di Base	805771740VARPE
Product Code	GB0459.AC222
Intended use	All clinical areas
Manufacturer	GARDHEN BALANCE SRL
Total weight	22 kg
Maximum load evenly distributed	50 kg



Sizes



Gardhen Bilance S.r.l.

Via G. Luraghi c/o Consorzio Il Sole - lotto S - 80038

Pomigliano d'Arco (NA) - Italia

Tel. +39 081 8692160 (selezione passante 5 linee)

E-mail: info@gardhenbilance.it

PEC: gardhenbilance@pec.it



Securtop *light*

Monobloc countertop

The images reproduced in this sheet are purely indicative and the dimensional values are nominal. It is the manufacturer's right to make changes to the designs without prior notice. For the identification of standard equipment and the selection of any additional components, please refer to the optional components catalogue and/or the sales proposals. In the event of discrepancies, the content of the sales proposals shall prevail.

For normative references please consult the user manual and declaration of conformity)

Dep. 437 rev. 18/02/2026