CHAIRSACCESSORIES




Accessory compatibility table

- 1 Stephen H Anatomical
- 2 Stephen H Anatomical Mini
- **3** Stephen H Anatomical Junior
- 4 Stephen H Anatomical Bariatrica
- 6 Stephen H Flat
- 7 Stephen H Flat Junior
- 8 Stephen H Flat Bariatrica
- **9** Stephen H Freedom
- 10 Stephen H Freedom Junior
- 11 Stephen H Freedom Bariatrica
- **12** Stephen H Comfort 2.0 3.0
- 13 Stephen H Comfort 2.0 3.0 Bariatrica
- 14 Stephen H Comfort 4.0
- **15** Stephen H Comfort 5.0
- **16** Stephen H Comfort 6.0
- 17 Stephen H Comfort 7.0
- 19 Stephen H Comfort 9.0
- 20 Stephen H Eco versione 1
- 21 Stephen H Eco versione 1 Junior
- 22 Stephen H Eco versione 2
- 23 Stephen H Eco versione 2 Junior
- 24 Stephen H Eco 2.0
- 25 Unica
- **26** Athena 2.0

- **27** Athena
- 28 Home Chair
- 29 Stephen H Ybrid
- **30** Stephen H Easy
- **31** Up
- 32 Akka
- **33** Stephen Scale H Anatomical
- **34** Stephen Scale H Anatomical Mini
- **35** Stephen Scale H Anatomical Junior
- **36** Stephen Scale H Anatomical Bariatrica
- 37 Stephen Scale H Flat
- **38** Stephen Scale H Flat Junior
- **39** Stephen Scale H Flat Bariatrica
- **40** Stephen Scale H Freedom
- 41 Stephen Scale H Freedom Junior
- **42** Stephen Scale H Freedom Bariatrica
- 43 Stephen Scale H Comfort 2.0 3.0
- 44 Stephen Scale H Comfort 2.0 3.0 Bariatrica
- **45** Athena Scale
- **46** CM/6R
- 47 Stephen Scale
- 48 Stephen H Infinity 1.0
- 49 Stephen


EXAMPLE OF ACCESSORY CONFIGURATION



Once the device has been configured, consult the Gardhen Bilance sales offices in order to avoid any incompatibilities between the chosen accessories.

The images reproduced are purely indicative and the manufacturer has the right to make changes to the projects without notice.

AC2P




Serving table on accessory rod equipped with height and rotation adjustment. Easily removable. Dimensions 234x185 mm.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	19	20	21	22	23	24	27	29
31	33	34	35	36	37	38	39	40	41	42	43	44	45	47	48	49							

AC3P



Hook for drainage bags also available on DIN rail.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	19	20	21	22	23	24	26	27
29	30	31	33	34	35	36	37	38	39	40	41	42	43	44	45	47	48	49					

AC4P



Removable, height-adjustable **adjustable IV pole** by mechanism.

Plastic hooks - available with 2 and 4 hooks

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	19	20	21	22	23	24	27	29
30	31	33	34	35	36	37	38	39	40	41	42	43	44	45	47	48	49						

AC4P/5



Removable, height-adjustable **adjustable IV pole** by mechanism.

Stainless steel hooks - available with 2 and 4 hooks

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	19	20	21	22	23	24	27	29
30	31	33	34	35	36	37	38	39	40	41	42	43	44	45	47	48	49						

AC4P/1



Removable, height-adjustable **non adjustable IV pole** by mechanism.

Plastic hooks - available with 2 and 4 hooks

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	19	20	21	22	23	24	27	29
31	33	34	35	36	37	38	39	40	41	42	43	44	45	47	48	49							

AC4P/6



Removable, height-adjustable **non adjustable IV pole** by mechanism.

Stainless steel hooks - available with 2 and 4 hooks

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	19	20	21	22	23	24	27	29
31	33	34	35	36	37	38	39	40	41	42	43	44	45	47	48	49							

AC4P/2




Height-adjustable **IV stand on wheels** with 4 hooks.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	19	20	21	22	23	24	25	26
27	28	29	30	31	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49		

AC4P/3



IV pole on DIN rail.


Incompatible

with the Stephen H Comfort 3.0 and Stephen Scale H Comfort 3.0..

COMPATIBLE FOR THE FOLLOWING CHAIRS

6 7 8 12 13 22 26 28 37 38 39 43 44

AC4P/4



Shaped steel IV pole with 4 hooks..

AC4P/S



Predisposition for IV pole placed on base (adjustable)..

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	19	20	21	22	23	24	26	27
28	29	30	31	33	34	35	36	37	38	39	40	41	42	43	44	45	47	48	49				

AC5P



Twin wheels diameter 100 mm.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	19	20	21	22	23	24	27	29
31	33	34	35	36	37	38	39	40	41	42	43	44	45	47	48	49							

AC5P/1



Twin wheels diameter 150 mm.

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	20	21	22	23	24	26	27	29
33	34	35	36	37	38	39	40	41	42	43	44	45	47	49									

AC6P



10-key pushbutton panel with spiral cable.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	20	21	22	23	24	26	27	29
33	34	35	36	37	38	39	40	41	42	43	44	45	47	49									

AC6P/2



Trendelenburg adjustment. Pedal controls located on the base.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	20	21	22	23	24	26	27	29
33	34	35	36	37	38	39	40	41	42	43	44	45	47	49									

AC6P/3



Synchronized backrest and leg adjustment.

Pedal controls located on the base. (Input/Exit Position)..

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	20	21	22	23	24	26	27	29
33	34	35	36	37	38	39	40	41	42	43	44	45	47	49									

AC6P/4



Seat inclination adjustment. Pedal controls located on the base. (Configurable with AC35 accessory).

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	20	21	22	23	24	26	27	29
33	34	35	36	37	38	39	40	41	42	43	44	45	47	49									

AC6P/5



Variable height adjustment. Pedal controls located on the base. (Configurable with AC28 accessory).

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	20	21	22	23	24	26	27	29
33	34	35	36	37	38	39	40	41	42	43	44	45	47	49									

AC6P/6



Backrest adjustment. Pedal controls located on the base.

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	20	21	22	23	24	26	27	29
33	34	35	36	37	38	39	40	41	42	43	44	45	47	49									

AC6MP



13-key control panel with 3 memory buttons, with spiral cable.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	20	21	22	23	24	26	27	33
34	35	36	37	38	39	40	41	42	43	44	45	47	49										

AC7P1G



Simultaneous brake on rear wheels in galvanized steel.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	20	21	22	23	24	26	27	29
33	34	35	36	37	38	39	40	41	42	43	44	45	47	49									

AC7P1S



Simultaneous brakeon rear stainless steel wheels...

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	20	21	22	23	24	26	27	29
33	34	35	36	37	38	39	40	41	42	43	44	45	47	49									

AC7P2G



Centralized brake on 4 wheels with directional system and activation via rear bilateral pedal, in galvanized steel.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	19	26	27	29	31	33	34	35
36	37	38	39	40	41	42	43	44	45	47	48	49											

AC7P2S




Centralized brake on 4 wheels with directional system and operation via rear bilateral pedal, in stainless steel.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	26	27	29	31	33	34	35	36
37	38	39	40	41	42	43	44	45	47	48	49												

AC7P3



Bilateral brake pedal front and rear.

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	19	26	27	29	31	33	34	35
36	37	38	39	40	41	42	43	44	45	47	48	49											

AC7P4



Centralized motorized brake, operated by pedal control.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1 4 6 8 9 11 12 13 14 15 16 17 27 29 33 36 37 39 40 42 43 44 45 47

49

AC7P5G



Centralized brake without directional wheels.

COMPATIBLE FOR THE FOLLOWING CHAIRS

	1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	26	27	29	33	34	35	36	37
3	8	39	40	41	42	43	44	45	47	49														

AC8P



Quick release of the backrest using the emergency lever (CPR).

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	19	20	21	22	23	26	27	31
33	34	35	36	37	38	39	40	41	42	43	44	45	47	48	49								

AC9B



Shaped and sewn armrest cover made of 50 g/m2 TNT

Pack of 2 pieces, minimum order quantity 100 pieces.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	19	20	21	22	23	24	25	26
27	28	29	30	31	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49		

AC9P



Shaped and sewn armchair cover made of 50 g/m2 TNT. Pack of 1 piece, minimum order quantity 60 pieces.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	19	20	21	22	23	24	25	26
27	28	29	30	31	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49		

AC10P



Shock position: backrest parallel to the floor, seat/legs inclined depending on the choice of the AC35 series..

COMPATIBLE FOR THE FOLLOWING CHAIRS

1 2 3 6 7 9 10 12 14 15 16 17 20 22 27 33 34 35 37 38 40 41 43 45 47 49

AC11



Roll holder positioned on the back of the armchair, in stainless steel for rolls with a diameter of 140/600 mm.

COMPATIBLE FOR THE FOLLOWING CHAIRS

6 7 8



Roll holder positioned on the back of the armchair, in stainless steel for rolls with a diameter of 140/600 mm.

With standard armrests.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1 3 4



Roll holder positioned on the back of the armchair, in stainless steel for rolls with a diameter of 140/600 mm

With multifunction armrests.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1 2 3 4



Rollholder ppositioned on the back of the armchair, in stainless steel for rolls with a diameter of 140/600 mm.

COMPATIBLE FOR THE FOLLOWING CHAIRS

9 10 11

AC11/1



Roll holder positioned on the back of the armchair, in stainless steel for rolls with a diameter of 140/600 mm.

After-sales supply.

List price does not include assembly.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1 3 4 6 7 8 9 10 11 20 21 22 23 24 27 33 35 36 37 38 39 40 41 42

45 47 49

AC11/C



Roll holder positioned on the back of the armchair, in stainless steel for rolls with a diameter of 140/600 mm.

COMPATIBLE FOR THE FOLLOWING CHAIRS

12 | 13 | 14 | 15 | 16 | 17 | 31 | 43 | 44 | 48



AC12




Independent movements of the backrest and leg sections. Additional engine to the standard equipment.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	19	20	21	22	23	24	25	26
27	28	29	30	31	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	49			

AC13



Motorized telescopic footrest with HPL shelf and support made of steel, adjustment via pushbutton panel.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	20	21	22	23	27	29	33	34	35	36	37	38	39	40
41	42	45	47	49																			

AC13/1



Motorized telescopic footrest equipped with further manual extension with HPL shelf and support made of steel, adjustment via pushbutton panel

1	2	3	4	6	7	8	9	10	11	20	21	22	23	27	29	33	34	35	36	37	38	39	40
4	42	45	47	49																			

AC13/C1







Poggiapiedi motorizzato in painted steel. The surface is reversible and therefore adaptable for patients of all heights. Adjustment via pushbutton. Easily removable.

Support surface with padding that can be created at different heights according to needs. The padding is covered with wear-resistant imitation leather material.

COMPATIBLE FOR THE FOLLOWING CHAIRS

12 | 13 | 14 | 43 | 44 |

AC13/C2








Motorized stainless steel footrest. The surface is reversible and therefore adaptable for patients of all heights. Adjustment via pushbutton. Easily removable.

Support surface with padding that can be created at different heights according to needs. The padding is covered with wear-resistant imitation leather material.

COMPATIBLE FOR THE FOLLOWING CHAIRS

12 | 13 | 14 | 43 | 44

AC13/SH1



Motorized footrest in steel and equipped with further manual extension

1	2	3	4	6	7	8	9	10	11	20	21	22	23	26	27	33	34	35	36	37	38	39	40
41	42	45	47	49																			

AC14P



Basket with stainless steel lid with revolving and height-adjustable support, positioned on the accessory rod. Dimensions 250x120x60h mm

On request it is available with different dimensions

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	19	20	21	22	23	24	27	29
31	33	34	35	36	37	38	39	40	41	42	43	44	45	47	48	49							

AC14P/1



Basket with removable and washable **stainless steel lid**, positioned on the armrest. Dimensions 250x120x60h mm

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	19	20	21	22	23	24	26	27
29	31	33	34	35	36	37	38	39	40	41	42	43	44	45	47	48	49						

AC15/PLUS





Serving table with double shelf made of ABS. The frame, on which the swivel wheels that allow movement are installed, is made of painted steel.

The shelf can be tilted into a lectern position and height-adjustable using pneumatic mechanisms. The load capacity is 10 kg.

Dimensions of the first shelf 420x200 mm.

Dimensions of the second shelf 420x550 mm.

Minimum height 830 mm.

Maximum height 1100 mm.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	4	6	8	9	11	12	13	14	15	16	17	19	20	22	24	25	26	27	28	29	30	31	33
36	37	39	40	42	43	44	45	46	47	48	49												

AC15/TOP





Serving table with shelf made of ABS.

The frame, on which the swivel wheels that allow movement are installed, is made of painted steel.

The shelf can be tilted into a lectern position and height-adjustable using pneumatic mechanisms. The load capacity is 10 kg.

Shelf dimensions 420x550 mm. Minimum height 830 mm. Maximum height 1100 mm.

1	4	6	8	9	11	12	13	14	15	16	17	19	20	22	24	25	26	27	28	29	30	31	33
36	37	39	40	42	43	44	45	46	47	48	49												

AC16



Telescopic footrest with HPL top, with manually adjustable stainless steel support.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	20	21	22	23	24	27	33	34	35	36	37	38	39	40
41	42	45	47	49																			

AC16/1



Padding covered in imitation leather for AC16 footrest, ideal for optimizing posture without shoes. The padding is fixed to the HPL panel using special knobs.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	12	13	14	20	21	22	23	24	26	27	29	33	34	35
36	37	38	39	40	41	42	43	44	45	47	49												

AC17P



Additional battery box useful for guaranteeing the functionality of the chair even in the absence of power from the mains.

This battery, in addition to guaranteeing functionality in emergency cases, ensures numerous operating cycles thanks to its high autonomy.

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	20	21	22	23	24	27	29	31
33	34	35	36	37	38	39	40	41	42	43	44	45	47	49									

AC17P/A



N. 2 batteries with battery charger (4.5 Ah).

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	15	20	21	22	23	24	26	27	29	30	31	33	34	35
36	37	38	39	40	41	42	45	47	49														

AC19



Armrest mechanism made of stainless materials.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1 2 3 4 6 7 8 9 10 11 20 21 22 23 24 33 34 35 36 37 38 39 40 41 42

AC21/HP2



Elliptical side walls made of HPL.

Not compatible

with external armrests (AC30H1).

COMPATIBLE FOR THE FOLLOWING CHAIRS

AC25P



Support for button panel housing integrated under the armrest of the chair..

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	19	20	21	22	23	24	26	27
29	30	31	33	34	35	36	37	38	39	40	41	42	43	44	45	47	49						

AC26

ROUTINE THAT CONTROLS WEIGHT LOSS DURING THE DIALYTIC SESSION AND GENERATES ALARMS. SOFTWARE FUNCTION FOR DEVICES EQUIPPED WITH SCALE..

Weight loss alarms.

COMPATIBLE FOR THE FOLLOWING CHAIRS

33 34 35 36 37 38 39 40 41 42 43 44 45 46 47

AC27



Additional headrest cushion adjustable in height, with padding made of expanded resin and antibacterial imitation leather covering.

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	19	20	21	22	23	24	26	27
31	33	34	35	36	37	38	39	40	41	42	43	44	45	47	48	49							

AC27/1



Additional headrest cushion adjustable in height, with padding made of expanded resin and antibacterial imitation leather covering.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	19	20	21	22	23	24	26	27
31	33	34	35	36	37	38	39	40	41	42	43	44	45	47	48	49							

AC27/2



Additional headrest cushion adjustable in height, with padding made of expanded resin and antibacterial imitation leather covering.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	19	20	21	22	23	24	26	27
31	33	34	35	36	37	38	39	40	41	42	43	44	45	47	48	49							

AC27/3



Additional headrest cushion radjustable in height, with padding made of expanded resin and antibacterial imitation leather covering.

MEDIUM DENSITY

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	19	20	21	22	23	24	26	27
31	33	34	35	36	37	38	39	40	41	42	43	44	45	47	48	49							

AC27/4



Additional headrest cushion adjustable in height, with padding made of expanded resin and antibacterial imitation leather covering.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	19	20	21	22	23	24	26	27
31	33	34	35	36	37	38	39	40	41	42	43	44	45	47	48	49							

AC27/5



Additional headrest cushion adjustable in height, with padding made of expanded resin and antibacterial imitation leather covering.

COMPATIBLE FOR THE FOLLOWING CHAIRS

29

AC27/6



Additional headrest cushion adjustable in height, with padding made of expanded resin and antibacterial imitation leather covering.

LOW DENSITY

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	19	20	21	22	23	24	26	27
31	33	34	35	36	37	38	39	40	41	42	43	44	45	47	48	49							

AC27/7



Additional headrest cushion adjustable in height, with padding made of expanded resin and antibacterial imitation leather covering.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	19	20	21	22	23	24	26	27
31	33	34	35	36	37	38	39	40	41	42	43	44	45	47	48	49							

AC27/8



Additional headrest cushion adjustable in height, with padding made of expanded resin and antibacterial imitation leather covering.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	19	20	21	22	23	24	26	27
31	33	34	35	36	37	38	39	40	41	42	43	44	45	47	48	49							

AC27/9



Additional headrest cushion adjustable in height, with padding made of expanded resin and antibacterial imitation leather covering.

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	19	20	21	22	23	24	26	27
31	33	34	35	36	37	38	39	40	41	42	43	44	45	47	48	49							

AC28/200



Variable height adjustment via motorized telescopic column with rectangular section with a <u>capacity of 200 kg and an excursion of 300 mm</u>.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1 2 3 6 7 9 10 12 13 14 15 16 17 26 27 33 34 35 37 38 40 41 43 44 45 47 49

AC28/200-248



Variable height adjustment via motorized telescopic column with rectangular section with a <u>capacity of 200 kg and an excursion of 248 mm</u>

COMPATIBLE FOR THE FOLLOWING CHAIRS

24

AC28/200-400



Variable height adjustment via motorized telescopic column with rectangular section with a <u>capacity of 200 kg and an excursion of 400 mm.</u>

COMPATIBLE FOR THE FOLLOWING CHAIRS

1 | 6 | 9 | 12 | 14 | 15 | 16 | 17 | 19 | 27 | 31 | 33 | 37 | 40 | 43 | 45 | 47 | 49

AC28/300



Variable height adjustment via motorized telescopic column with rectangular section with a capacity of 300 kg and an excursion of 250 mm. For bariatric chairs.

COMPATIBLE FOR THE FOLLOWING CHAIRS

4 8 11 13 29 36 39 42 44

AC29



Armrest band made of imitation leather, adjustable via Velcro.

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	19	20	21	22	23	24	26	27
28	29	30	31	33	34	35	36	37	38	39	40	41	42	43	44	45	48						

AC30P



Multifunction armrestin stainless steel and aluminum adjustable in: height, adjustable, tilting, folding and removable.



COMPATIBLE FOR THE FOLLOWING CHAIRS

1 2 3 4 6 7 8 9 10 11 20 21 22 23 27 33 34 35 36 37 38 39 40 41 42 45

AC30/C1














Ergonomic multifunction armrest, adjustable via knob.

Movements allowed:

- integral with the backrest(foto 1)
- swivel (foto 2)
- adjustable in height (foto 3)
- removable (foto 4)
- foldable (foto 5)
- tiltable "up/down" (foto 6).

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	19	20	21	22	23	24	27	31
33	34	35	36	37	38	39	40	41	42	43	44	45	48										

AC30/C2



Implementation of the AC30/C1 armrest with the **function tion slide**.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	19	20	21	22	23	24	27	31
33	34	35	36	37	38	39	40	41	42	43	44	45	48										

AC30/C3



Implementation of the AC30/C1 armrest with the **internal rotation function**.

	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	19	20	21	22	23	24	27	31
3	3 34	35	36	37	38	39	40	41	42	43	44	45	48										

AC30/C5



Slim multifunction armrest in steel adjustable via knob with slide function.

Armrest dimensions: length 515 mm, width 105 mm..

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	19	20	21	22	23	24	27	31
33	34	35	36	37	38	39	40	41	42	43	44	45	48										

AC30/H1



Multifunctional armrest adjustable in height, swivel and removable if necessary.

COMPATIBLE FOR THE FOLLOWING CHAIRS

29

AC31



Joint cushion rmade with expanded resin and antibacterial imitation leather covering.

Slim multifunction armrest in steel adjustable via knob with slide function.

Thanks to its ergonomics and its cylindrical shape, it is functional to optimize comfort during long stays..

COMPATIBLE FOR THE FOLLOWING CHAIRS

1 2 3 4 6 7 8 20 21 22 23 27 33 34 35 36 37 38 39 4!





Courtesy **LED light** equipped with "ON/OFF" easily accessible to the patient

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	20	21	22	23	27	29	33	34
35	36	37	38	39	40	41	42	43	44	45	47	49											

AC33



Manually **height-adjustable headrest** with 4 intermediate positions. Mechanism made of stainless material and easily removable if necessary.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1 2 3 4 6 7 8 19 20 21 22 23 27 31 33 34 35 36 37 38 39 45 48

AC33/1



Reclining headrest 0/15°.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1 | 4 | 6 | 8 | 33 | 36 | 37 | 39 | 47 | 49





Headrest reclining 0/15° and **adjustable in alsize** manually.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1 | 4 | 6 | 8 | 33 | 36 | 37 | 39 | 47 | 49

AC34P



Key for inhibiting push button controls. This device, supplied to healthcare personnel, can be used to protect patient safety.

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	19	20	21	22	23	24	25	26
27	28	29	30	31	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	49			

AC35



Fixed seat parallel to the floor.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1 2 3 4 33 34 35 36

AC35/1



Motorized seat with adjustable inclination from 0° to -15°, operated via pushbutton. This function is useful for passive postural variation of the patient for the prevention of bedsores.

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	19	20	22	24	26	27	29	31
33	34	35	36	37	38	39	40	41	42	43	44	45	47	49									

AC35/2



Motorized seat with adjustable inclination from +3° to -12°, operated via pushbutton. This function is useful for passive postural variation of the patient for the prevention of bedsores...

COMPATIBLE FOR THE FOLLOWING CHAIRS

1 2 3 6 7 9 10 20 22 24 27 33 34 35 37 38 40 41 45 47 49

AC35/3



Motorized seat with adjustable inclination from -10° to -25°, operated via pushbutton. This function is useful for passive postural variation of the patient for the prevention of bedsores.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1 | 2 | 3 | 6 | 7 | 9 | 10 | 20 | 22 | 24 | 27 | 33 | 34 | 35 | 37 | 38 | 40 | 41 | 45 | 47 | 4

AC35/4



Motorized seat with adjustable inclination from -12° to +12°, operated via pushbutton.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1 2 3 33 34 35 47 49

AC35/5



Motorized seat with adjustable inclination from 0° to +10°, operated via buttons.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1 2 3 33 34 35 47 49



Pair of transparent PVC covers to protect the armrests, useful for preserving the integrity of the upholstery over time.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	19	20	21	22	23	24	27	28
29	30	31	33	34	35	36	37	38	39	40	41	42	43	44	45	47	48	49					

AC36/1



Transparent PVC covers for footrests, useful for preserving the integrity of the upholstery over time.

1	2	3	4	6	7	8	9	10	11	12	13	14	20	21	22	23	24	27	28	29	30	33	34
35	36	37	38	39	40	41	42	43	44	45	47	49											

AC36/2



Transparent PVC cover for the leg section integral with the cushion, useful for preserving the integrity of the upholstery over time.

On request it can also be made removable.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1 2 3 4 6 7 8 19 20 21 22 23 26 27 28 30 31 33 34 35 36 37 38 39

45 47 48 49

AC37



Extra large armrest 200 mm.

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	20	21	22	23	24	27	33	34
35	36	37	38	39	40	41	42	43	44	45													

AC37/1



Bracciolo extra long 625x140 mm.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1 2 3 4 6 7 8 9 10 11 20 21 22 23 27 33 34 35 36 37 38 39 40 41

42 45

AC37/3



Elongated armrest made with expanded resin on the back.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1 2 3 4 6 7 8 9 10 11 20 22 27 33 34 35 36 37 38 39 40 41 42 45

AC37/4



Seamless armrest supplement.

COMPATIBLE FOR THE FOLLOWING CHAIRS

 1
 2
 3
 4
 6
 7
 8
 9
 10
 11
 12
 13
 14
 15
 16
 17
 20
 21
 22
 23
 27
 33
 34
 35

 36
 37
 38
 39
 40
 41
 42
 43
 44
 45

AC37/A



Extra large armrest 360x120 mm.

COMPATIBLE FOR THE FOLLOWING CHAIRS

32

AC38





Additional cushion for lowered seat:

- Padding with standard height (photo 1).
- Seat in the lowered version (photo 2).

COMPATIBLE FOR THE FOLLOWING CHAIRS

1 | 2 | 3 | 4 | 6 | 7 | 8 | 20 | 21 | 22 | 23 | 27 | 33 | 34 | 35 | 36 | 37 | 38 | 39 | 4



Leg section at 90° to the floor.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1 | 2 | 3 | 4 | 6 | 7 | 8 | 27 | 33 | 34 | 35 | 36 | 37 | 38 | 39 | 45

AC40



Air-adjustable internal lumbar inflatable cushion. Functional to optimize posture during the phases of therapy with preventive antidecubitus effect.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	12	13	14	20	21	22	23	27	33	34	35	36	37	38
39	40	41	42	43	44	45	47	49															

AC40/1



Air-adjustable external lumbar inflatable cushion. Functional to optimize posture during the phases of therapy with preventive antidecubitus effect.

1	2	3	4	6	7	8	9	10	11	12	13	14	20	21	22	23	24	26	27	28	29	30	33
34	35	36	37	38	39	40	41	42	43	44	45	47	49										

AC40/3



Additional external lumbar cushion. Functional to optimize posture during the phases of therapy with a preventive anti-decubitus effect.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	12	13	14	20	21	22	23	24	26	27	28	29	30	33
34	35	36	37	38	39	40	41	42	43	44	45	47	49										

AC41



Thermal printer for armchairs with weighing system, useful for printing the weight value in real time.

It is powered by rechargeable batteries and the supplied power supply, and can also be fixed to the chair with a dedicated support...

COMPATIBLE FOR THE FOLLOWING CHAIRS

33 34 35 36 37 38 39 40 41 42 43 44 45 46 47

AC43P



Sick lift pole made of reinforced steel tubing and ergonomic triangular handle, adjustable in height via a special mechanism.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1 | 4 | 6 | 8 | 9 | 11 | 20 | 22 | 27 | 29 | 33 | 36 | 37 | 39 | 40 | 42 | 45 | 47 | 49



TV 24" installed on an armchair using a dedicated swivel support.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	12	13	14	20	21	22	23	27	33	34	35	36	37	38
39	40	41	42	43	44	45	47	49															

AC48/PP



Preparation for TV support..

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	12	13	14	20	21	22	23	27	33	34	35	36	37	38
39	40	41	42	43	44	45	47	49															

AC48/SP



TV stand.

										12	13	14	20	21	22	23	27	33	34	35	36	37	38
1	39	40	41	42	43	44	45	47	49														




Pair of side panels made of HPL, functional to preserving the integrity of the base as well as performing its aesthetic function.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1 | 2 | 3 | 6 | 7 | 9 | 10 | 16 | 17 | 27 | 33 | 34 | 35 | 37 | 38 | 40 | 41 | 45

AC50P



Cable winding support

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	20	21	22	23	24	27	29	30
31	33	34	35	36	37	38	39	40	41	42	43	44	45	47	49								

AC51



Backrest made of memory foam.

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	20	21	22	23	27	33	34	35
3	37	38	39	40	41	42	43	44	45	47	49												



Seat made of memory foam

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	20	21	22	23	24	27	33	34
35	36	37	38	39	40	41	42	43	44	45	47	49											

AC55P



Built-in servant's table removable if necessary. Dimension 420x550 mm.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1 4 6 8 9 11 12 13 14 20 22 27 33 36 37 39 40 42 43 44 45

AC56



Handle made of stainless steel to facilitate transport activities.


COMPATIBLE FOR THE FOLLOWING CHAIRS

1 4 6 8 9 11 20 22 24 27 28 29 30 33 36 37 39 40 42 45 47 49





Handles made of PVC to facilitate transport activities.



COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	20	21	22	23	27	28	30	33
34	35	36	37	38	39	40	41	42	43	44	45	47	49										

AC56/C



Handle made of stainless steel to facilitate transport activities.

COMPATIBLE FOR THE FOLLOWING CHAIRS

12 | 13 | 14 | 16 | 17 | 19 | 43 | 44 | 48

AC56/C1



Handle made of aluminum to facilitate transport activities.

COMPATIBLE FOR THE FOLLOWING CHAIRS

12 | 13 | 14 | 15 | 16 | 17 | 31 | 43 | 44



Print for AC49 side panels to protect the base.



COMPATIBLE FOR THE FOLLOWING CHAIRS

1 2 3 6 7 9 10 27 33 34 35 37 38 40 41 45

AC60



Adjustable belts to protect the patient's safety

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	20	21	22	23	27	29	30	33
34	35	36	37	38	39	40	41	42	43	44	45	47	49										



Leg section at 90° to the seat

COMPATIBLE FOR THE FOLLOWING CHAIRS

1 2 3 6 7 16 17 19 20 22 27 28 31 33 34 35 37 38 45 48

AC63/1




Leg section at 90° to the seat for bariatric chairs.

COMPATIBLE FOR THE FOLLOWING CHAIRS

4 8 36 39

AC74



Automatic emergency button to obtain <u>Trendelenburg or shock positions</u>.

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	20	22	24	27	29	33	34	35
36	37	38	39	40	41	42	43	44	45	47	49												

AC74/1



Button to obtain automatic input position.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	19	20	22	27	31	33	34	35
36	37	38	39	40	41	42	43	44	45	47	49												

AC74/2



Button to obtain the clinostati position.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	20	22	27	29	33	34	35	36
37	38	39	40	41	42	43	44	45	47	49													

AC74/3



Emergency stop button.

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	19	20	22	26	27	29	31	33
34	35	36	37	38	39	40	41	42	43	44	45	47	49										



Adjustable cup holder located under the armrest.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	19	20	21	22	23	24	27	31
33	34	35	36	37	38	39	40	41	42	43	44	45	48										

AC75/1



Adjustable courtesy table placed under the armrest with slots for glasses or smartphones.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1 4 6 8 11 12 13 19 20 22 24 27 31 33 36 37 39 42 43 44 45 48

AC75/2




Adjustable cup holder placed under the armrest, extra-long version..

COMPATIBLE FOR THE FOLLOWING CHAIRS

1 | 4 | 6 | 8 | 11 | 12 | 13 | 19 | 20 | 22 | 24 | 27 | 31 | 33 | 36 | 37 | 39 | 42 | 43 | 44 | 45 | 48

AC75/3



Multifunctional adjustable support with hooks located under the armrest, useful for holding small bags.

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	20	21	22	23	24	27	31	33
34	35	36	37	38	39	40	41	42	43	44	45	48											



Additional cushion for back section.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	20	21	22	23	24	27	33	34
35	36	37	38	39	40	41	42	43	44	45	47	49											

AC79



Built-in thermo-cushion for backrest or seat.


1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	20	21	22	23	27	33	34	35
36	37	38	39	40	41	42	43	44	45	47	49												

AC81P



Articulated and removable servant's table, iin stainless steel and aluminum with storage edge. Made of non-porous material without bacterial retention.

Top dimensions equal to 440x340 mm.



COMPATIBLE FOR THE FOLLOWING CHAIRS

1 | 4 | 6 | 8 | 9 | 11 | 12 | 14 | 20 | 22 | 27 | 33 | 36 | 37 | 39 | 40 | 42 | 43 | 44 | 45 | 47 | 49


AC81/SP



Arrangement for serving table accommodation.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1 | 4 | 6 | 8 | 9 | 11 | 12 | 14 | 20 | 22 | 27 | 33 | 36 | 37 | 39 | 40 | 42 | 43 | 44 | 45 | 47 | 49



Leg section with motorized telescopic extension, useful in the case of configuration with 90° leg section to give greater comfort to the patient.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1 6 27 33 37 45

AC82/1



Leg section with motorized telescopic extension for bariatric chairs.

COMPATIBLE FOR THE FOLLOWING CHAIRS

4 8 36 39

AC83




Cover made of ABS for base. Functional to protect and preserve the sensitive parts of the device.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1 | 6 | 9 | 12 | 13 | 14 | 15 | 26 | 27 | 33 | 37 | 40 | 43 | 45 | 47 | 49

AC83/C



Cover made of ABS for base. Functional to protect and preserve the sensitive parts of the device.

COMPATIBLE FOR THE FOLLOWING CHAIRS

16 | 17 | 19 | 31 | 48

AC83/H



Cover made of ABS for base. Functional to protect and preserve the sensitive parts of the device.

COMPATIBLE FOR THE FOLLOWING CHAIRS

4 8 11 13 29 44

AC84



Cover made of ABS for seat. Functional to protect and preserve the sensitive parts of the device.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1 6 9 20 22 27 33 37 40 45 47 49

AC84/C



Cover made of ABS for seat section. Functional to protect and preserve the sensitive parts of the device.

COMPATIBLE FOR THE FOLLOWING CHAIRS

12 14 15 16 17 19 31 43 48

AC84/H



Cover made of ABS for seat section. Functional to protect and preserve the sensitive parts of the device.

COMPATIBLE FOR THE FOLLOWING CHAIRS

29

AC85/C



Cover made of ABS for leg section. Functional to protect and preserve the sensitive parts of the device.

COMPATIBLE FOR THE FOLLOWING CHAIRS

12 | 14 | 15 | 43

AC86/C



Cover made of ABS for backrest section.

Functional to protect and preserve the sensitive parts of the device.

COMPATIBLE FOR THE FOLLOWING CHAIRS

12 14 15 16 17 19 31 43

AC86/H



Cover made of ABS for backrest section.

Functional to protect and preserve the sensitive parts of the device.

COMPATIBLE FOR THE FOLLOWING CHAIRS

29


AC87



Bi-elastic silicone coating, extremely resistant and stain-resistant. To choose the color, consult the appropriate color card.

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	19	20	21	22	23	24	25	26
27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	

AC87/1



Leatherette covering with stain-resistant properties. To choose the color, consult the appropriate color card.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	19	20	21	22	23	24	25	26
27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	

AC88



Armrest cover made of material: waterproof, washable, breathable, bi-elastic and bacteriostatic.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	19	20	21	22	23	24	27	28
29	30	31	33	34	35	36	37	38	39	40	41	42	43	44	45	48							

AC89



Armchair cover made of material: waterproof, washable, breathable, bi-elastic and bacteriostatic.

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	19	20	21	22	23	24	26	27
28	29	30	31	33	34	35	36	37	38	39	40	41	42	43	44	45	47	48	49				

AC90/H



Armchair cover in TNT material.

COMPATIBLE FOR THE FOLLOWING CHAIRS

29

AC90/H1



Cotton armchair cover.

COMPATIBLE FOR THE FOLLOWING CHAIRS

29

AC94P



Wireless call kit <u>installed on the armrest of the chair</u>.

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	20	21	22	23	24	27	28	29
30	33	34	35	36	37	38	39	40	41	42	43	44	45										

AC95P




Cylinder holder installed on the base of the chair.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1 4 6 8 9 11 12 13 14 15 16 17 20 22 26 27 33 36 37 39 40 42 43 44

45 47 49

AC96P




Wireless call button installed on the armrest of the chair.

COMPATIBLE FOR THE FOLLOWING CHAIRS

26

AC99



90° backrest section.

COMPATIBLE FOR THE FOLLOWING CHAIRS

 1
 2
 3
 4
 6
 7
 8
 9
 10
 11
 12
 13
 14
 15
 16
 17
 20
 21
 22
 23
 27
 33
 34
 35

 36
 37
 38
 39
 40
 41
 42
 43
 44
 45
 47
 49



Adjustable side rails installed on DIN bar.

Incompatible

with the Stephen H Comfort 3.0 and Stephen Scale H Comfort 3.0.

COMPATIBLE FOR THE FOLLOWING CHAIRS

12 13 26 28 37 39 43 44

AC109/1



Adjustable side rails installed on DIN bar.

Incompatible

with the Stephen H Comfort 3.0 and Stephen Scale H Comfort 3.0.



Leg rests installed on a DIN bar adjustable using a special mechanism.

COMPATIBLE FOR THE FOLLOWING CHAIRS

6 8 22 37 39

AC115



Activation of the controls for moving the sections on the Fly Plus display:

backrest, legs, variable height and memory positions: entrance position and emergency position.

COMPATIBLE FOR THE FOLLOWING CHAIRS

33 34 35 36 37 38 39 40 41 42 43 44 45 46 47

AC116



Flexible document holder positioned on the back of the document chair. Size 260x350 mm.

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	19	20	21	22	23	24	25	26
27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	



Multiple transport handle for Athe-na 2.0, useful for transport both in a postural position and in a stretcher position. These handles facilitate maneuverability of the device even with the patient on board.

COMPATIBLE FOR THE FOLLOWING CHAIRS

26

AC118



5ª directional wheel with pedal activation to facilitate movement for long journeys from place to place.

COMPATIBLE FOR THE FOLLOWING CHAIRS

26

AC119



Pelvic girdle.

COMPATIBLE FOR THE FOLLOWING CHAIRS

12 | 13 | 14 | 15 | 16 | 17 | 26 | 43 | 44



Independent pneumatic movement of the backrest and legrest sections.

COMPATIBLE FOR THE FOLLOWING CHAIRS

32

AC132



Compass sides with servo-assisted folding mechanism made of polymer and aluminium.

COMPATIBLE FOR THE FOLLOWING CHAIRS

26

AC133



Seat width extendable by 20 cm on DIN rail.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1 4 6 8 9 11 12 13 14 15 16 20 22 26 27 29 33 36 37 39 40 42 43 44

45



Bilateral steel DIN bar for housing additional accessories measuring 25x8 mm.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	4	6	8	9	11	12	13	14	15	16	20	24	22	26	27	29	33	36	37	39	40	42	43
44	45	47	49																				

AC135



Support feet with non-slip base.

COMPATIBLE FOR THE FOLLOWING CHAIRS

32

AC136



N. 2 rear directional wheels diameter 125 mm equipped with anti-slip front tips.



Multifunction armrest with radial joint installed on DIN bar.

COMPATIBLE FOR THE FOLLOWING CHAIRS

26

AC145



Support with radial joint on DIN bar and IV pole.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1 4 6 8 9 11 12 13 14 15 16 20 22 26 27 29 33 36 37 39 40 42 43 44 45 47 49

AC146



USB port to allow charging/powering of electronic devices with safety cover.

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	20	21	22	23	24	25	26	27
28	29	30	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	49				



Leg section equipped with footrest with folding HPL top.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1 | 4 | 6 | 8 | 20 | 22 | 27 | 33 | 36 | 37 | 39 | 45

AC153



Footrest with swivel mechanism and automatic braking system.

COMPATIBLE FOR THE FOLLOWING CHAIRS

16 | 17


AC157



Headrest with 3D adjustment mechanism integrated into the back. Eight different types of neck supports available.

COMPATIBLE FOR THE FOLLOWING CHAIRS

15 | 17 | 3



Column stand for battery chargers.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1	2	3	4	6	7	8	9	10	11	12	13	14	15	16	17	19	20	21	22	23	24	25	26
27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	

AC164



90° rotatable display for scale chairs.

COMPATIBLE FOR THE FOLLOWING CHAIRS

33 34 35 36 37 38 39 40 41 42 43 44 45 46 47

AC167



GBTouch multifunction control touch screen display.

COMPATIBLE FOR THE FOLLOWING CHAIRS

33 34 35 36 37 38 39 40 41 42 43 44 45 46 47





WI-FI/Bluetooth combo module for GBTouch.

COMPATIBLE FOR THE FOLLOWING CHAIRS

33 34 35 36 37 38 39 40 41 42 43 44 45 46 47

AC167/2




GBTouch **USB** stick for storage

COMPATIBLE FOR THE FOLLOWING CHAIRS

33 34 35 36 37 38 39 40 41 42 43 44 45 46 47

AC167/3



10/100 ethernet port for GBTouch.

COMPATIBLE FOR THE FOLLOWING CHAIRS

33 34 35 36 37 38 39 40 41 42 43 44 45 46 47

AC167/4



USB A port for GBTouch. Weight data transmissions. Not compatible with the AC167/2 USB stick.

COMPATIBLE FOR THE FOLLOWING CHAIRS

33 34 35 36 37 38 39 40 41 42 43 44 45 46 47

AC167/5




USBBport for GBTouch. External printer connection.

COMPATIBLE FOR THE FOLLOWING CHAIRS

33 34 35 36 37 38 39 40 41 42 43 44 45 46 47

AC167/6




RS-232 port for GBTouch.

COMPATIBLE FOR THE FOLLOWING CHAIRS

33 34 35 36 37 38 39 40 41 42 43 44 45 46 47

AC167/7



RS-485 port for GBTouch.

COMPATIBLE FOR THE FOLLOWING CHAIRS

33 34 35 36 37 38 39 40 41 42 43 44 45 46 47

AC174



Legrest with leg section 90°/180° with respect to the floor and with retractable folding footrest.

COMPATIBLE FOR THE FOLLOWING CHAIRS

1 6 33 37

Once the device has been configured, consult the Gardhen Bilance sales offices in order to avoid any incompatibilities between the chosen accessories.

The images reproduced are purely indicative and the manufacturer has the right to make changes to the projects without notice.



Gardhen Bilance S.r.l.

Via G. Luraghi c/o Consorzio II Sole lotto S - 80038 Pomigliano d'Arco (NA) - Italy Tel. +39.081.8692160 pbx Fax +39.081.8692460 E-mail: info@gardhenbilance.it PEC: gardhenbilance@pec.it