

Home Chair

The right chair when it comes to narrow spaces



CHAIR

Frequent usage

Chemotherapy

Infusions therapy

Blood donor

Blood samples

Emergency

Dialysis

Multifunctional





10 BUTTONS HANDSET
WITH 2 MEMORY FUNCTION

Benefits

- Anti-tilt frame, it prevents any tilt in every operative conditions.
- Frame painted with epoxy paints enriched with antibacterial. Available in a wide range of colours.
- Ergonomic and extremely comfortable surface.
- Padding in non-deformable polyurethane foam.
- Two separate cushions for back and seat/leg. The cushions are designed to be ergonomic and they help to prevent pressure ulcers.
- The chair can be equipped with a wide range of accessories (please refer to the accessories catalogue).
- Leatherette upholstery with the following features: extremely resistant, antibacterial, flame retardant class 11M, antimycotic, It is also resistant to: alcohol, hydrogen peroxide, hypochlorite sodium max 5% and to commonly used disinfectants.
- Max load capacity 180 kg.

Standard features

- Painted steel frame in thermosetting powder enriched with antibacterial. Available in a wide range of colours.
- Lying surface divided in 3 sections: backrest, seat and leg.
- Simultaneous and independent movements of back and leg sections through 2 motors - Equipment code AC12.
- 90° leg section compared to the seat - Equipment code AC63.
- No. 4 wheels Ø 100 mm in stainless steel (2 rear wheels with single brake).
- Chair's movements are controlled by handset with a 3,5mtl coiled cable. The handset is user friendly and has visual pictograms useful for partially sighted people. Handset is equipped by 10 buttons with 2 memories. For patient's safety purpose the handset is provided with an inhibition key - Equipment code AC34P. Insulation IP66.
- Independent and removable cushions.
- Upholstery in flame retardant leatherette UNI 9175/87 e 9175FA-1/94 class 11M - Equipment code AC87/1.




Movements

OVERTURNING ARMREST WITH AUTOMATIC ROTATIVE DEVICE THAT FOLLOWS THE INCLINATION OF THE BACKREST ENSURING A COMFORTABLE STAY TO THE PATIENT.

90° LEG SECTION COMPARED TO THE SEAT



TRENDELENBURG -10°



MOTORIZED AND SYNCHRONIZED EXTENSION / ELEVATION OF THE LEG SECTION

















Colour chart





























FRAME COLOURS





















The entire structure of the chair is protected by epoxy and thermosetting paints. The paints are enriched with antibacterial upon request. Whether not specified the frame of the device will be painted in white 9010.

 Code 2004	 Code 5012	 Code 7035
 Code 3020	 Code 6027	 Code 9005
 Code 1018	 Code 6018	 Code 9010
 Code 1021	 Code 7001	 Code 4010
 Code 5010	 Code 9006	

COLOURS FOR STANDARD UPHOLSTERY



 Code 5947	 Code 0872	 Code 0014
 Code 5948	 Code 4947	 Code 0027
 Code 0680	 Code 1944	 Code 0006
 Code 0508	 Code 4008	 Code 0346
 Code 0041	 Code 1946	 Code 016043
 Code 0064	 Code 4933	 Code 011499
 Code 0055	 Code 0430	 Code 015815
 Code 1816	 Code 0151	 Code 015814
 Code 0833	 Code 0667	 Code 016042

COLOURS FOR SILICONIC UPHOLSTERY (THEY HAVE A HIGHER PRICE) ACCESSORY CODE AC87

 Codice 7039.08	 Codice 7039.14	 Codice TH.439
 Codice 7039.06	 Codice TH.513	 Codice TH.438
 Codice 7039.17	 Codice TH.511	 Codice TH.322
 Codice 7039.20	 Codice TH.305	 Codice TH.224
 Codice 7039.10	 Codice TH.378	 Codice 312
 Codice 7039.19	 Codice TH.300	 Codice TH.226
 Codice 7039.13	 Codice TH.510	



Dimensions



MAIN CORD



Type F - Schuko



Type L - 16A



Type I - Australia

Type B - Canada
and USA

Type G - UK plug

Technical sheet



GMDN Code	38447
ID Code	GB0138.SP Painted steel GB0138.SS Stainless steel
Intended use	Multifunctional chair
Manufacturer	GARDHEN BILANCE SRL
Overall weight	70 kg
Wheels	No. 4 100 mm stainless and swivable (rear wheels are equipped with brakes)
Power	100 - 240 Vac 50/60 Hz
Electric motors	Low voltage (24 V)
Max electrical input	4A
Insulation	IP44
Max load	180 kg
Safety working load	215 kg

For certifications please refer to "Certifications and Quality system" section

Pictures contained in this brochure is for illustrative purpose only. The manufacturer reserves the right to modify the projects without notice. For accessories identification please refer to the accessory catalogue or commercial offer. When data are different please refer to commercial offer.

Dep. EN_84 rev. 07/12/2021

