

Stephen H Flat – Junior

Essential design



CHAIR

Frequent usage

- Chemotherapy
- Infusions Therapy
- Day Hospital
- Blood donor
- Blood samples
- Emergency
- Dialysis
- Patient's transport
- Dermatology / Cosmetics
- Multifunctional
- Reduced footprint
- Paediatrics





ARMRESTS MOVEMENT IS SYNCHRONIZED WITH THE BACKREST

Benefits

- Structure with a reduced footprint.
- Anti-tilt frame, it prevents any tilt in every operative conditions.
- Frame painted with epoxy paints enriched with antibacterial. Available in a wide range of colours.
- Simultaneous movement of back and leg section through 1 motor; it can be implemented with 2 up to 5 motors according to the requested configuration:
 - Dedicated motor for the independent movement of the backrest and the legrest (accessory code AC12).
 - Motorized column for height adjustment (accessory code AC28).
 - Dedicated motor for the inclination of the seat (accessory codes AC35/1 - AC35/2 - AC35/3).
 - Dedicated motor for footrest adjustment (accessory code AC13 - AC13/1); (please refer to the commercial offer and the accessory catalogue;
- Three separate cushions for back, seat and leg functional for a perfect ergonomomy and posture preventing decubitus wounds.
- The upholstery is completely removable allowing a thorough and effective sanitization reducing the risks of cross infections. A wide range of colours is available.
- Upholstery's features: antibacterial, flame retardant class 11M, antimycotic, ecologic, latex free, phtalate free and resistant.
- Armrest's movements synchronized with the backrest.
- Excellent and ergonomic posture with decubitus wound preventive effect.
- Braking system ensures high stability to the device.
- Easily sanitizable with the most common used detergents.
- Armrest with multiple adjustment functions.
- Max load capacity 150 kg.



THE CUSHIONS ARE FULLY REMOVABLE ALLOWING AN EASY SANITIZATION OF THE DEVICE. SURFACE IN NON-POROUS HPL PANELS.

Standard features

- Painted steel frame in thermosetting powder. Wide range of colours at page 5.
- Lying surface divided in 4 sections (headrest, backrest, seat and leg).
- Simultaneous movement of the back section and the leg section through 1 motor.
- 10° Trendelenburg position (shock position as alternative).
- Headrest joint to the backrest.
- Armrest's movements synchronized with the backrest. The armrests are completely removable and can be rotate, tilted and overturned at 90°; (please refer to the accessories catalogue for other version of armrests).
- No. 4 wheels Ø 100mm antistatic, sturdy, stainless and swivable.
- Simultaneous brake on rear wheels in chromed steel - Equipment code AC7P1G.
- Chair's movements are controlled by handset with a 3,5mtl coiled cable. The handset is user friendly and has visual pictograms useful for partially sighted people. Handset is equipped by 10 buttons with 2 memories. For patient's safety purpose the handset is provided with an inhibition key - Equipment code AC34P. Insulation IP66.
- Bilateral handset support - Equipment code AC25P.
- Side panels in HPL - equipment code AC49.
- Easily removable soft foam independent cushions.
- Upholstery in flame retardant leatherette UNI 9175/87 e 9175FA-1/94 class 11M - Equipment code AC87/1



**10 BUTTONS HANDSET
WITH 2 MEMORIES**

Armrests' movements



ROTATION



TILTABLE




OVERTURNABLE




JOINTLY TO THE BACKREST

Movements

TRENDELENBURG



SHOCK POSITION



RECLINER POSITION



Colour chart



FRAME COLOURS

The entire structure of the chair is protected by epoxy and thermosetting paints. The paints are enriched with antibacterial upon request. Wether not specified the frame of the device will be painted in white 9010.

Code 2004	Code 5012	Code 7035
Code 3020	Code 6027	Code 9005
Code 1018	Code 6018	Code 9010
Code 1021	Code 7001	Code 4010
Code 5010	Code 9006	

COLOURS FOR STANDARD UPHOLSTERY



Code 5947	Code 0872	Code 0014
Code 5948	Code 4947	Code 0027
Code 0680	Code 1944	Code 0006
Code 0508	Code 4008	Code 0346
Code 0041	Code 1946	Code 016043
Code 0064	Code 4933	Code 011499
Code 0055	Code 0430	Code 015815
Code 1816	Code 0151	Code 015814
Code 0833	Code 0667	Code 016042

COLOURS FOR SILICONIC UPHOLSTERY (THEY HAVE A HIGHER PRICE) ACCESSORY CODE AC87

Code 7039.08	Code 7039.14	Code TH.439
Code 7039.06	Code TH.513	Code TH.438
Code 7039.17	Code TH.511	Code TH.322
Code 7039.20	Code TH.305	Code TH.224
Code 7039.10	Code TH.378	Code 312
Code 7039.19	Code TH.300	Code TH.226
Code 7039.13	Code TH.510	



Dimensions



Technical sheet

GMDN Code	38447
ID Code	GB0002.JP Painted steel GB0002.JS Stainless
Inteded use	Multifunctional chair
Manufacturer	GARDHEN BILANCE SRL
Overall weight	70 kg
Wheels	No.4 Ø 100 mm swivable and stainless with braking system
Power	100 - 240 Vac 50/60 Hz
Electric motors	Low voltage 24 (Volt)
Max electrical input	4A
Insulation	IP44
Max load	150 kg
Safety working load	185 kg
Backrest angle	76° (customizable upon request)
Leg angle	45° (customizable upon request)

MAIN CORD



Type F - Schuko



Type L - 16A



Type I - Australia



Type B - Canada and USA



Type G - UK plug

For certifications please refer to "Certifications and Quality system" section

Pictures contained in this brochure is for illustrative purpose only. The manufacturer reserves the right to modify the projects without notice. For accessories identification please refer to the accessory catalogue or commercial offer. When data are different please refer to commercial offer.

Dep. EN_213 rev. 06/12/2021